


**Operating Manual** 

# Kangaroo<sup>™</sup>

ePump<sup>™</sup> Enteral Feed and Flush Pump with Pole Clamp, Programmable



## en

Thank you for purchasing the Kangaroo<sup>™</sup> ePump<sup>™</sup> Enteral Feed and Flush Pump with Pole Clamp. With proper care, this device will provide you with years of precision service.

## **Table of Contents**

P	a	a	e
	-	_	-

Section I — General Information	1
Section II — Safety and Warnings	3

Section III — Icon Identification	

## Section IV — Initial Setup

Attaching the A/C Power Adapter	1
Battery Setup7	'
Attaching Pole Clamp	'

## Section V — Instructions for Use

Quick Start	
General Startup	
Placement/Mounting	
A/C Power Operation	
Battery Power Operation	
Power On/Off	
Language Selection, First Power Up	
Keep or Clear Prior Pump Settings	
Loading Pump Sets	
Prime Pump	
Auto Priming	
Hold-To-Prime	
Feed-Only Pump Sets11	
Feed & Flush Pump Sets12	
Re-priming after Bag Empty12	
Selecting Feed Mode12	
EZMODE Option12	
Selecting EZMODE	

## **Table of Contents**

Selecting Continuous or Intermittent Feed Mode
Continuous Mode (Adjust Feed)
Continuous Mode Feeding14
Continuous Mode Flushing14
Intermittent Mode (Adjust Bolus)
Intermittent Mode Feeding15
Intermittent Mode Feeding - Bolus Max15
Intermittent Mode Flushing
Run
Keep Tube Open (KTO)
Lock Screen Feature
Hold17
Clear Volume Fed
Adjust Settings
Resume In Min
Run Mode
View History
Flush Now
More Options
Buzzer Alarm Volume
History
Language Selection
Continuous/Intermittent Mode
BioTech Menu
Accessing the BioTech Screen
BioTech Screen Data
BioTech Options
Lock Settings Feature (Program Locked)20
EZMode
Communication

## Section VI — Performance Evaluation

System Performance Tests	. 21
Certifying the Accuracy of Pump Flow Rate	. 21

Section VII — Cleaning		. 22
------------------------	--	------

## **Table of Contents**

Section VIII —	- Battery Replacement .	
----------------	-------------------------	--

## Section IX — Error/Warning/Informational Screens and Troubleshooting

System Error	26
Hold Error	28
Rotor Error	
Feed Error	28
Flush Error	29
Flow Error	29
Pump Set Dislodged Error.	
Battery Low	
Feeding Complete	
Pump Set Use >24 Hours Warning	
LED Array	

## Section X — Specifications and Symbols

Specifications	
Symbols on Pump and Accessories	

Section XI — Customer Service	34
-------------------------------	----

## Section XII — Maintenance

Side Door Enclosing Pump Set Loading Area	
Battery Pack	
Battery Door	
A/C Power Adapter	
Pole Clamp	
Section XIII — Service Part Numbers	35
Section XIV — Warranty	
Section XV — Electromagnetic Conformity Declaration	

This product contains software solely owned by Covidien. Covidien grants the user a non-exclusive, limited license to use the software according to the operating instructions. A copy of the license can be obtained from Covidien.

## Section I – General Information

The Kangaroo<sup>M</sup> ePump<sup>M</sup> enteral feed and flush pump with pole clamp is a simple-to-use, precision enteral feeding pump. It is easily programmed to provide patients with either continuous or intermittent feeding, and can provide automatic flushing capability when used with Kangaroo<sup>M</sup> ePump<sup>M</sup> feed & flush sets.

#### Intended use:

Intended for use in patients with any condition requiring enteral feeding and/or enteral hydration, which can be accomplished by means of an enteral feeding pump and pump set. The pump and feeding sets are intended to be used in alternate and acute care settings by clinicians. The purpose of this device is to deliver enteral nutrition at a controlled rate to a patient's gastrointestinal system.

#### **User Interface:**

- Intuitive user interface
- Large, backlit LCD display
- Step-by-step prompts to guide programming and operation
- "Stoplight" LED array visually indicates pump status in a bright or darkened room

#### **Ergonomics**:

- Quiet operation
- Compact, aesthetic design
- Tabletop usage or IV pole mounting
- · Simple loading of pump sets
- Transparent side door to protect pump set
- Integrated handle for easy transport

#### Features:

- MISTIC (Magnetic Intelligent Set Type Identification Connector) Feed-Only vs. Feed & Flush Set identification system to ensure match between the pump's user interface and pump set type
- Automatic Anti-Free-Flow (AFF) System
- · Audible alarm to indicate errors or pump set loading conditions
- · Attitude-independent design eliminates need for drip chamber on pump set
- · Sensor technology detects both upstream and downstream flow conditions
- · Continuous feed, intermittent feed, and flush capability
- Auto-prime feature reduces the need for time-consuming manual priming
- Hold feature with restart and Keep Tube Open (KTO) feature
- Stores previous 72 hours of feeding and flushing history

#### Service:

- Serial I/O port for data retrieval
- Uses replaceable, rechargeable NiMH Battery Pack for 15 hours of backup power
- Removable A/C Power Adapter

Kangaroo<sup>™</sup> ePump<sup>™</sup> Enteral Feed and Flush Pump with Pole Clamp




Figure 1A. Kangaroo™ ePump™ features, front and rear views.




Figure 1B. Kangaroo™ ePump™ pump set loading area.

## Section II — Safety and Warnings

**Caution:** U.S. federal law restricts the sale of this device to physicians or to their direct representatives.

- 1. Read this booklet thoroughly before using the Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump.
- 2. Do not use this device near flammable anesthetics.
- 3. Disconnect power supply before cleaning or servicing.
- 4. Use only Kangaroo<sup>™</sup> ePump<sup>™</sup> Enteral Feed and Flush Pump A/C Power Adapter with pump. See Section XIII Service Part Numbers for replacement of power adapter. Use of an alternate consumer style A/C power adapter may cause damage to the charging circuit and battery of the Kangaroo<sup>™</sup> ePump<sup>™</sup> feeding pump.
- 5. **Danger:** Use only Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed or feed with flush pump sets with this device. Pump is not compatible with other pump sets. Use of other feeding sets with this pump can create hazardous situations, including free-flow conditions that can result in overfeeding, underfeeding, formula in the lungs, and death to a patient.
- 6. For performance tests, see Section VI Performance Tests. For other integrity checks, consult with a qualified biomedical technician or contact the manufacturer (Section XI Customer Service).
- 7. For service or for technical information, please contact Customer Service (Section XI).
- 8. Do not open the main housing, as there are no user-serviceable parts inside. Opening may affect function of device and voids the warranty. To replace battery, see Section VIII.
- 9. Dispose of old battery-powered electronic equipment in a manner consistent with institutional policy for expired equipment disposal.
- 10. Cleaning frequency and practices must be consistent with institutional policy for cleaning of non-sterile devices. See Section VII Cleaning, for instructions on cleaning the Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump.

- 11. See icon descriptions in Section III Icon Identification and Section X Specifications and Symbols for additional safety information.
- 12. This device is designed to minimize the effects of uncontrolled electromagnetic interference and other types of interference from external sources. Avoid use of other equipment that may cause erratic operation or degradation in the performance of this device.
- 13. **Warning:** The serial I/O port (see Figure 1A) is intended only for use with the approved communication cable (Part Number 382499).

Any other equipment that can be connected to the DB9 Serial I/O port may lead to a reduced level of safety of the pump.

- The Covidien approved Communication Cable (Part Number 382499) is specifically approved as complying with the safety requirements for accessories used with the pump, whereas the safety approval of the accessory has been performed in accordance to the appropriate IEC 60601-1 and/or IEC 60601-1-1 harmonized standard.
- Connection of the pump to other equipment could result in previously unidentified risks. Before connecting to other equipment, the risks of using the pump with this equipment should be identified, analyzed, evaluated, and controlled prior to use.
- 14. **Caution:** This enteral feeding pump should only be used for patients who can tolerate the flow rates and accuracy levels delivered by the pump. Premature infants may require higher accuracy rates than specified for this enteral feeding pump.
- 15. Do not attempt to rotate the valve stem inside the feeding set valve assembly. Rotation of the valve stem can only be done by the pump while the valve assembly is properly loaded into the pump.
- 16. Should feeding sets require rinsing, it is recommended that the feeding sets be rinsed while they are loaded in the pump.
- 17. Used feeding sets should be disposed of in accordance with current hospital procedure or local disposal guidelines.
- 18. **Caution:** Ensure buzzer hole is unobstructed during normal operation so as to allow clear recognition of alarm. Inability to hear the alarms could pose a serious risk to the patient, since the operator may not hear an alarm.
- 19. **Caution:** The battery cells used in this device may present a fire or chemical hazard if mistreated. Do not disassemble, heat above 100°C (212°F), or incinerate.
- 20. **Danger:** Strangulation Hazard. Avoid leaving power adapter wires, feeding set tubing or other choking hazards where infants or young children can become caught. If these objects get wrapped around a child's neck, strangulation and death can occur.
- 21. **Danger:** The pump and disposable feeding set all contain small parts which could become detached and pose a choking hazard. Some of these components could be inhaled or swallowed by a small child, toddler, or infant, which could result in suffocation and death. Keep all small components out of reach of small children.
- 22. **Danger:** Explosion Hazard. Do not use the pump in the presence of flammable anesthetics. Flammable anesthetics can ignite due to a spark within the unit, which could result in fire or explosion.
- 23. **Danger:** Do not power the pump while plugged into a power outlet above 2000m altitude, especially while on an airplane or other air vehicle. Do not attempt to charge the pump battery above 2000m altitude. Failure to comply could result in overheating, fire, or explosion of the internal rechargeable battery pack. Run the pump only using the internal pump battery power in these situations.
- 24. **Caution:** Do not store the pump or power adapter at temperatures >50°C (122°F). This can damage the equipment sensors, which will prevent the pump from operating under normal conditions.
- 25. **Warning:** Do not modify this equipment without authorization of the manufacturer. Modification of any devices or accessories can result in physical hazards including delayed therapy, over delivery, under delivery, electrocution, and fire. These hazards could result in patient injury or death.
- 26. **Caution:** The power adapter wires, feeding set tubing, and pump accessories may cause a tripping hazard. Avoid leaving wires, cords, or tubing in a pathway where a person could trip on them and sustain an injury.

- 27. **Caution:** Avoid using accessories, detachable parts and materials with the pump that are not recommended in this manual. Use only approved Kangaroo<sup>™</sup> ePump<sup>™</sup> accessories with the pump.
- 28. **Caution:** Use the pump only as directed in this user manual. Do not interconnect this device with other devices or modify the equipment in any way outside of the recommendations in this manual. Failure to comply could result in incorrect delivery of formula to the patient and could result in damage to the pump.
- 29. **Caution:** This pump is not intended to be used in MRI environments or in the presence of strong magnetic fields. Do not use these devices in any areas with strong magnetic fields. The pump contains metal components which could cause unintended movement. This unexpected movement could cause harm due to falling objects or collisions.
- 30. **Caution:** There are significant hazards associated with accidental misconnections with other infusion devices, which could lead to patient harm or death. For more information about hazards and risk reduction strategies associated with misconnections, see the following: The Joint Commission Sentinel Event Alert Issue 36 April 13, 2006.
- 31. Do not use the pump for delivery of any fluids or substances that are not enteral solutions prescribed by qualified medical personnel.
- 32. The pump is designed to be used outdoors for short periods of time (no more than 24 hours). Leaving the pump outdoors for extended periods of time (exceeding 24 hours) can result in damage and/or fading of the pump devices.
- 33. This device is designed for use on a conventional IV pole. As with any medical device, it is possible for the weight of the pump to cause the IV pole to tip over. This could result in injury to a patient or operator. When attaching the pump to the IV pole, take precautions to ensure the IV pole remains stable while in use.
- 34. This enteral feeding system was designed to meet IEC 60601-1 safety standards. For clarification, purposes, the feeding set is considered an Applied Part and has been tested and evaluated accordingly.
- 35. **Warning:** Not for intravenous use. Do not use for intravenous infusion into a patient. Intravenous infusion of enteral fluids can result in serious complications up to and including death.
- 36. **Warning:** This enteral feeding pump should only be used for patients who can tolerate the flow rates and accuracy levels delivered by the pump. Premature infants may require higher accuracy rates than specified for this enteral feeding pump. Delivery of fluid to patients who cannot tolerate the pump accuracy can result in over or under delivery with the possibility of aspiration.
- 37. Use only commercially available pre-packed or commercially prepared feeding solutions prescribed by a licensed health care provider, dietitian or nutritionist. Do not use homemade blenderized or liquidised foods or other non-prescribed, non-commercially available feeding solutions.

#### **Run Mode**

Vertically-moving droplet in *RUNNING* screen indicates Run mode.

## **Battery Symbols**

- Full charge
- 2/3 or less charge
- 1/3 or less charge
- Very low battery (15 minutes max.)

## **Auto Priming Indicator**

When auto priming, the bar will progress to show auto prime feature is active.

## **Volume Setting**

Indicates alarm volume setting.

## **Error Indicator**

Pump operations are stopped until the error is resolved.

## **Information Indicator**

This indicates information only, and does not require any immediate action.

## EZMODE Indicator

Indicates that EZMODE feature is turned on.

## **Communication Indicator**

Indicates that communication feature is turned on.

## Volume to be Delivered Indicator

Indicates that VTBD feature is turned on.

## Keep Tube Open

Indicates that Keep Tube Open feature is turned on.





## PRIMING


















## Section IV — Initial Setup

## Attaching the A/C Power Adapter

Plug the A/C Power Adapter (Part Number 382491) into the A/C adapter port on the back of the pump. See Figure 1A to locate port on pump. Make note of the location of the pins and tab on the plug to properly orient the plug.

The connector, which attaches to the back of the pump, has a spring-loaded locking sleeve that locks into the pump, thus helping to prevent the connector from being accidentally unplugged from the pump. To remove the connector, grasp the outer sliding sleeve member and pull. Do not detach the connector by pulling on the cord only.



adapter port.

## **Battery Setup**

The battery pack is fully charged before shipping, and is left unconnected to optimally preserve battery life.

**Warning:** The battery should be connected before using the pump.

It is also recommended that the battery charge be "topped off" before battery-only operation is utilized. The battery pack will charge continuously whenever the pump is plugged into a wall outlet, including during normal use of the pump with A/C power. Eight hours of charging is required to fully recharge the battery pack.

A new battery pack will provide approximately 15 hours of service before needing to be recharged.

## **Attaching Pole Clamp**

The Kangaroo<sup>TM</sup> ePump<sup>TM</sup> enteral feed and flush pump can be attached to a vertical pole using the pole clamp, included with the pump. The pole clamp provides the ability to retain the power cord to prevent accidental cord removal or cord loss. Simply place the A/C adapter power cord into the U-shaped groove located on the back of the pole clamp as shown in Figure 3 and attach the pole clamp to the Kangaroo<sup>TM</sup> ePump<sup>TM</sup> enteral feed and flush pump. Be sure not to pinch the cord of the A/C power adapter or damage the cord insulation when installing.

To attach the pole clamp to the Kangaroo<sup>TM</sup> ePump<sup>TM</sup> enteral feed and flush pump, align the hole on the pole clamp with the mounting hole on the back of the pump. Use a 3/8'' - 16 bolt to fasten together, as shown in Figure 3.

If cord restraint is preferred but use of the pole clamp is not desired, a cord restraint clip is provided with the Kangaroo<sup>M</sup> ePump<sup>M</sup> enteral feed and flush pump. Using the pole clamp mounting bolt, follow the pole clamp assembly instruction to attach the cord clip (Refer to Fig. 3).

Figure 3. Attaching and using the pole clamp or cord clip.

## Section V — Instructions For Use

### **Quick Start**

- 1. Press **U POWER** on lower, right-hand corner of control panel.
- 2. Open blue door enclosing pump set loading area.
- 3. To load pump set , refer to Figure 4:
  - Grasp finger tab on valve and insert into left pocket (1).
  - Grasp black ring retainer and gently stretch tubing around rotor (2). Pull up on retainer and insert into right pocket (3).
  - Press down on finger tab to ensure that valve is fully seated. Finger tab should align with locating rib.
  - Close blue door.

(Note: display status line should read SET LOADED).

- 4. To automatically prime the pump, press ► **PRIME PUMP** and ► **AUTO PRIME**. For Kangaroo<sup>™</sup> ePump<sup>™</sup> feed & flush sets, auto priming will prime both lines.
- 5. To directly control pump priming, press ► **PRIME PUMP** and then press and hold the **HOLD TO PRIME** button(s) ► . If a feed & flush set is used, remember to first press ► **HOLD TO PRIME FLUSH** until the flushing fluid has reached the valve, and then press ► **HOLD TO PRIME FEED** until the feeding fluid has reached past the valve and down to the distal end connector.
- 6. Use buttons ▶ on left side of LCD to set the feeding parameters. Set flushing parameters if a feed & flush set is loaded. The continuous or intermittent feed mode selection is made under the **MORE** options menu.
- 7. Press > RUN when ready to start. The screen will display RUNNING.
- 8. To stop, press **HOLD** or press and hold the **OPOWER** button.




Figure 4. Installation of pump set.

## **General Startup**

## **Placement/Mounting**

There are two recommended placement methods for the Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump.

- Attached to a vertical IV pole via the pole clamp peripheral device, included with the pump.
- Placed feet down on any horizontal surface.

## **A/C Power Operation**

Plug the Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump into an A/C power outlet for normal operation. When the pump is not plugged in, or if A/C power is interrupted, a built-in rechargeable battery will run the pump.

## **Battery Power Operation**

If A/C power is not available to the pump, or if A/C power is interrupted, the pump will automatically operate on backup battery power. The built-in battery recharges any time the pump is plugged into an outlet.

A fully-charged, new battery will supply 15 hours of backup power at 125 mL/hr before requiring recharging. Contact Customer Service, Section XI, to replace the battery pack if battery life degrades.

#### Power On/Off

To power up the pump, press the **D POWER** button on the lower right of the front panel. To power down the pump, press and hold the **D POWER** button. A countdown will begin on the display; when the countdown reaches 0, immediately release the POWER button. **Note:** The display backlight and Status LED array will remain on for a short period after releasing the POWER button while the pump completes its internal power down process and will then extinguish.

## Language Selection, First Power Up

Immediately after powering up the pump for the first time, the language option screen will appear. Select the language using the ( $\uparrow$ ) or ( $\checkmark$ ) buttons. When the language is selected, push the fifth button ( $\rightarrow$ ) to enter the language selection. After the language is selected the first time the pump is powered up, the language selection screen will not appear during future power up operations. Should the language need to be changed after the initial power up, it can be changed through the "**MORE OPTIONS**" menu (see Figure 10A).

#### **Keep or Clear Prior Pump Settings**

Immediately after powering up the pump, if the prior feed or flush settings were not cleared and the settings are not locked in the BioTech Screen, the pump will give the following two options:

**KEEP SETTINGS** - Select this option to start with the same settings that were most recently programmed into the pump. Any of the previously-programmed settings can be modified, if needed. The feed volume totals, "**mL FED**" (Continuous Feed Mode) or "**BOLUSES FED**" (Intermittent Feed Mode) and "**mL FLUSH**" are also maintained.

**CLEAR SETTINGS** - Select this option to reset all input feeding settings and the output feed totals to zero. It will then be necessary to program all settings before running the pump. The BioTech options settings will not change.

The status line of the pump will then show *LOAD A SET* (Figure 5), or, if the pump set is already loaded, it will show *SET LOADED* (Figure 7).

**Note:** The output feed totals ("**mL FED**," "**BOLUSES FED**," "**mL FLUSH**") can be cleared at any time, and without clearing the input settings, by doing the following:

- Start the pump **RUNNING**
- Press 
  HOLD
- Press 
  CLEAR VOLUME

## **Loading Pump Sets**

The Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump will indicate **LOAD A SET** in the screen's status line if a pump set is not installed. A blinking icon of a pump set will also appear in the upper left corner of the screen. Figure 5A shows the screen for continuous feeding mode and Figure 5B shows the screen for intermittent feeding mode (intermittent feeding mode not available in EZMODE).

To load a pump set, follow the steps below:

- 1. Open the blue transparent door that encloses the rotor and pump set loading area.
- 2. Load the pump set by holding the finger tab on the valve, and inserting the valve down into the slot on the left. (Figure 6A)
- 3. Grasp the MISTIC retainer end (black ring retainer) and wrap the tubing counterclockwise around the rotor. **Avoid overstretching the silicone tubing. (Figure 6B)**
- 4. Carefully pull the MISTIC retainer end upwards to position it over the slot on the right, and then lower the retainer into the slot. (Figure 6B)
- 5. Press down on the valve's finger tab to ensure proper seating.
- 6. Close blue transparent door.
- 7. The bottom of the feed bag should be 18 inches from the top of the pump.




Figure 5A. *LOAD A SET* opening menu screen for Continuous Feeding mode (EZMODE not activated).




Figure 5B. *LOAD A SET* opening menu screen for Intermittent Feeding mode (EZMODE not activated).






Figure 6A.





Kangaroo<sup>™</sup> ePump<sup>™</sup> Enteral Feed and Flush Pump with Pole Clamp

The display should read **SET LOADED** and will show one of the following opening menus, depending on the type of pump set detected, and whether continuous or intermittent feeding mode has been selected.






Figure 7A. Opening menu for Continuous Feed mode with a feed-only set.

Figure 7B. Opening menu for Continuous Feed mode with a feed & flush set.



Figure 7C. Opening menu for Intermittent Feed mode with a feed-only set.





Figure 7D. Opening menu for Intermittent Feed mode with a feed & flush set.

## **Prime Pump**

The Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump may be primed automatically with a single button push, including the flushing line if a feed & flush set is loaded. The pump can also be primed in a more interactive method using the hold-to-prime feature. The flow rate for priming, whether Auto Prime or Hold-to-Prime, is 1960 ml/hr (32.7 mL/minute). **Note:** During priming the Flow Error alarm is disabled.

After a pump set has been loaded and the status line shows **SET LOADED**, press > PRIME PUMP in the opening menu (Figures 7A-7D) to get one of the **PRIME PUMP** menus shown in Figures 8A-8B.



feed-only set.




Figure 8B. PRIME PUMP menu for feed & flush set.

## **Auto Priming**

The **AUTO PRIME** option will be available if the pump senses that the currently installed pump set was not previously primed, and that there is no fluid in the line near the rotor. If auto-priming is desired, press and release > **AUTO PRIME** to automatically prime the line(s). For feed & flush pump sets, the auto prime feature will automatically prime both lines, starting with the flush line.

#### Press > **STOP** to cancel **AUTO PRIME**.

When auto-priming is completed, the status line will show **AUTO PRIME COMPLETE** and will no longer show the **AUTO PRIME** option. Be sure that the line(s) are fully primed. If not, use the hold-to-prime option(s) to complete the priming, as described below.

## Hold-to-Prime

The hold-to-prime options allow for precise interactive control of the priming process.

## **Feed-Only Pump Sets**

For feed-only pump sets, the screen will appear as in Figure 8A. Press and hold > HOLD TO PRIME FEED until the feed line has been primed down to the stepped end connector at the end of the pump set.

Kangaroo™ ePump™ Enteral Feed and Flush Pump with Pole Clamp

## Feed & Flush Pump Sets

For feed & flush pump sets, the display will appear as in Figure 8B. First press and hold > HOLD TO PRIME FLUSH to prime the flushing line, and then press and hold > HOLD TO PRIME FEED to prime the feed line down to the stepped end connector at the end of the pump set.

If a feed & flush pump set is used, remember to first press > HOLD TO PRIME FLUSH until the flushing fluid has reached the valve, and then press > HOLD TO PRIME FEED until the feeding fluid has reached past the valve and down to the stepped end connector. If the feeding line is primed first, subsequent priming of the flushing line will force the air between the flush solution bag and valve into the main line, which already has been primed with enteral formula.

## **Re-priming after Bag Empty**

A pump set bag that has been emptied will trigger the **FEED ERROR** screen. In this condition the pump set bag can be refilled to continue the feeding, but only after the pump set has been re-primed.

To re-prime the pump set, do the following:

- Disconnect the feeding line from the patient
- Refill the bag
- Press CONTINUE to begin the pump RUNNING
- Press > HOLD
- Press ADJUST SETTINGS
- Press 
  PRIME PUMP
- Press HOLD TO PRIME to prime the line(s)

Pump sets should not be reused after 24 hours of initial usage. The feeding set should also be replaced after 24 hours from initiation of feeding. This ensures that the system is operating within specified parameters and prevents bacterial growth that could be a hazard to the patient.

## **Selecting Feed Mode**

#### **EZMODE Option**

The **EZMODE** option is a limited, **CONTINUOUS MODE** only setting of ePump<sup>™</sup>. It feeds non-stop (no VTBD option) at a rate 0-400mL/hr as defined by the user until manually stopped. It has no **BOLUS**, **FLUSHING**, **RUN MODE SCREEN LOCK**, or **AUTORESUME** capability and displays no **HISTORY** information.

Turning on **EZMODE** will automatically force the unit out of the **INTERMITTENT (BOLUS FEED)** mode if set. Although the Continuous Mode feed VTBD will be set to 0 while **EZMODE** is enabled, the previous value of the **Feed VTBD** will be stored in flash memory and recalled when "**EZMODE**" is deselected. The user can change the Feed Rate by normal means. All other user settings, including those not compatible with the **EZMODE** (**BOLUS VTBD**, **NUMBER OF BOLUSES**, **SUPER BOLUS**, etc), that were set before the pump was put into the **EZMODE** will be ignored and invisible while in **EZMODE** but remain stored in the flash memory for when **EZMODE** is deselected. **EZMODE** essentially has no effect on **CLEAR SETTINGS**; pressing **CLEAR SETTINGS** while **EZMODE** is enabled will clear all feed and flush settings and totalizers as normal, not just those that apply to **EZMODE**. Selecting the **CLEAR SETTINGS** option on power up will not clear **Biotech Options** such as **EZMODE** or **LOCK SETTINGS**.


If the user selects the **ADJUST FEED** option from either the **LOAD A SET** or **SET LOADED** screens, control will pass directly to and from the Enter Feed Rate screen and bypass the **ADJUST FEED** screen since without the ability to adjust VTBD in **EZMode**, the Adjust Feed screen is redundant.

## Selecting EZMODE

**EZMODE** may be selected through the **Biotech Options** button which may be accessed on the Biotech screen (see page 19). The **EZMODE** setting will be stored in the flash memory and saved when the unit is shut off.

**Note:** If the user inserts a feed/flush set while the **EZMODE** is enabled, **EZMODE** will be cancelled and the unit will immediately revert back to the normal mode. If **EZMODE** is desired, the user must then re-enable **EZMODE** through the **Biotech Options** screen.

When enabled and not in **CERTIFICATION MODE**, a small **EZ** icon will be visible on the bottom of the screen near the battery icon. However, it will not appear in the **POWER DOWN**, **SYSTEM ERROR**, or **LOW BATTERY SCREENS**.



## **Selecting Continuous or Intermittent Feed Mode**

The Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump can be programmed to feed continuously or intermittently.

The continuous feeding mode will deliver the enteral nutrition at a steady rate, either until the programmed volume has been delivered or until the supply has been exhausted. The Intermittent feeding mode delivers boluses of enteral nutrition at programmed time intervals. The bolus volume and feed rate are also programmed.

To select continuous feed mode or intermittent feed mode, press > **MORE** on the opening menu, shown in Figures 5A-5B or Figures 7A-7D, to access the **MORE OPTIONS** menu, Figure 10A.

The **MORE** option is available either before a pump set has been loaded (Figures 5A-5B) or after the pump set is loaded (Figures 7A-7D). See subsection "More Options" for information on the other options in the **MORE OPTIONS** menu.

Select > CONTINUOUS/INTERMITTENT under the *MORE OPTIONS* menu, Figure 10A. Then press > CONTINUOUS MODE or > INTERMITTENT MODE from the SELECT MODE menu, Figure 10B. Press > DONE to exit.

LANGUAGES
CONTINUOUS/INTERMITTENT
Figure 10A. MORE OPTIONS menu.




Figure 10B. *SELECT MODE* menu for continuous or intermittent mode selection.

## Continuous Mode — Adjust Feed


## **Continuous Mode Feeding**

For **CONTINUOUS** feeding mode, the **ADJUST FEED** option is available in the opening menus (Figure 7A for feed-only sets or Figure 7B for feed & flush sets). Select this option to program the feeding parameters from the **ADJUST FEED** screen, Figure 11A.

Press **FEED RATE** in Figure 11A to define the rate of delivery. Press the **buttons** in Figure 11B to program the rate from 1 to 400 in increments of 1 mL/hr. Press **FETER** (top button) to exit the menu.

Press > **FEED VTBD** in Figure 11A to define the Volume To Be Delivered. Press the > buttons in Figure 11C to program the volume from 1 to 3000 mL in increments of 1 mL. Press > **ENTER** (top button) to exit the menu. If the **FEED VTBD** is set to zero, the pump will run until the supply is exhausted.

**Note:** When a **FEED VTBD** is set, the remaining volume to be delivered (remaining VTBD) will be displayed next to the v icon on the lower right-hand portion of the screen (Figure 11A). If the **FEED VTBD** is reprogrammed during a feeding, from **HOLDING** mode, the pump will restart the feeding and deliver the entire volume of the new VTBD setting.






Figure 11B. Setting the Feed Rate.




Figure 11C. Setting the Feed Volume to be Delivered.

## **Continuous Mode Flushing** (Not Available in EZMODE)

For **CONTINUOUS** feeding mode with a feed & flush set loaded, the opening menu will appear as shown in Figure 7B. Press > **ADJUST FLUSH** in Figure 7B to program the flushing parameters from the **ADJUST FLUSH** screen, shown in Figure 12A.


Press > **FLUSH VOLUME** in Figure 12A to define the volume per flush occurrence. Press the > buttons in Figure 12B to program the flush volume from 10 to 500 in increments of 1 mL. Press > **ENTER** (top button) to exit the menu.


Press **FLUSH INTERVAL** in Figure 12A to define the time interval between starts of flushing. Press the **buttons in Figure 12C to** program the time interval from 1 to 24 hours in increments of 1 hour. Press **ENTER** (top button) to exit the menu.

The pump will automatically limit flushing settings which exceed pump capabilities to deliver.

**Caution:** The Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump flushes at a rate of 1960 mL/hr (32.7 mL/minute). Use care when programming the flush volume so that it matches the patient's need.

**Note:** The pump will automatically flush an additional 25 mL of water after completing the programmed delivery of formula, as long as at least one flush had occurred during formula delivery.








Figure 12B. Setting the Flush Volume.




Figure 12C. Setting the Flush Time Interval.

Kangaroo™ ePump™ Enteral Feed and Flush Pump with Pole Clamp

## **Intermittent Mode** — Adjust Bolus (Not Available in EZMODE)

### **Intermittent Mode Feeding**

For **INTERMITTENT** feeding mode, the **ADJUST BOLUS** option is available in the opening menus (Figure 7C for a feed-only set or Figure 7D for a feed & flush set). Select this option to program the feeding parameters from the **ADJUST BOLUS** screen, Figure 13A.

Press > **BOLUS RATE** in Figure 13A to define the rate of delivery.

Press the ▶ buttons in Figure 13B to program the delivery rate from 1 to 400 in increments of 1 mL/hr. Press ▶ ENTER (top button) to exit the menu.

Press > **NUMBER OF BOLUSES** in Figure 13A to define the total number of boluses to be delivered.

Press the > buttons in Figure 13C to program the number of boluses, up to a maximum of 99.

Press > **ENTER** (top button) to exit the menu. The number of boluses can be set as zero, in which case the pump will run until the supply is exhausted.

Press > VOLUME PER BOLUS in Figure 13A to define the volume to be delivered in each bolus.

Press the ▶ buttons in Figure 13D to program the volume per bolus from 1 to 3000 mL in increments of 1 mL.

Press **• ENTER** (top button) to exit the menu.

Press **BOLUS INTERVAL** in Figure 13A to define the time interval between starts of bolus deliveries.

Press the ▶ buttons in Figure 13E to program the time interval between starts of bolus deliveries from 1 to 24 hours in increments of 1 hour. Press ▶ **ENTER** (top button) to exit the menu.

The pump will automatically limit bolus settings which exceed pump delivery capabilities. For example, a bolus delivery volume of 200 mL cannot be made once per hour if a rate of 150 mL/ hr is already programmed. In this case, the maximum volume allowed would be 150 mL.

**Note:** If the NUMBER OF BOLUSES is reprogrammed during a feeding (from HOLD mode), the pump will restart the bolus feeding and will deliver all the boluses as defined by the new setting.

## Intermittent Mode Feeding – Bolus Max (Not Available in EZMODE)

The **BOLUS MAX** option can be used to deliver a bolus at a high rate, similar to gravity feeding.

Press > BOLUS MAX in the BOLUS RATE menu, Figure 13B, to program the BOLUS MAX mode.

The feeding rate for **BOLUS MAX** is 999 mL/hr.




Figure 13A. ADJUST BOLUS menu.




Figure 13B. Setting the Bolus Rate.




Figure 13C. Setting the Number of Boluses.




Figure 13D. Setting the volume per Bolus.





Figure 13E. Setting the time interval between starts of bolus.

## **Intermittent Mode Flushing** (Feed flush set only)

For **INTERMITTENT** feeding mode with a feed & flush set, the opening menu, shown in Figure 7D, will show the **ADJUST FLUSH** option. Press > **ADJUST FLUSH** in Figure 7D to program the flushing parameters from the **ADJUST FLUSH** screen, Figure 14A.



**Note:** The nature of two intermittent processes, periodic bolus feeding and periodic flushing, is such that occasional conflicts can arise where each process needs to be running at the same time. The flushing process is assigned precedence over the feeding process for feed & flush pump sets. "Flushing precedence" means that a flush cycle will run to its end even if a bolus delivery is scheduled to begin. It also means that a flush cycle will interrupt a bolus feeding cycle that is already in progress. The feeding cycle will pick up where it was interrupted after the flushing cycle is complete. It is therefore possible that the net feeding volume is reduced from what may be intended, over a given time period.

Press > **FLUSH VOLUME** in Figure 14A to define the volume per flush cycle. Press the > buttons in Figure 14B to program the flush volume from 10 to 500 in increments of 1 mL. Press > **ENTER** (top button) to exit the menu.

Press > **FLUSH INTERVAL** in Figure 14A to define the time interval between starts of flushing. Press the > buttons in Figure 14C to program the time interval from 1 to 24 hours in increments of 1 hour. Press > **ENTER** to exit the menu.

The pump will automatically limit flushing settings which exceed pump capabilities to deliver.

**Caution:** The Kangaroo<sup>™</sup> ePump<sup>™</sup> feed and flush pump flushes at a rate of 1960 mL/hr (32.7 mL/minute). Use care when programming the flush volume so that it matches the patient's need.

**Note:** The pump will automatically flush an additional 25 mL of water after completing the programmed delivery of formula, as long as at least one flush had occurred during formula delivery.

## Run

The Kangaroo<sup>™</sup> ePump<sup>™</sup> Enteral feed and flush pump will only run when a pump set is properly loaded and the necessary feeding parameters have been programmed. For a feed & flush set, the flushing parameters do not have to be set for the pump to begin operating. If flushing parameters are set as zeros, there will be no flushing.

After the pump has been programmed, press > **RUN** to start operation. Figure 15A shows the **RUNNING** screen for continuous feeding with flushing. The status line will indicate **RUNNING**, and the droplet icon will move vertically on the screen. The green LED indicator will be illuminated to give a quick visual indication of positive status in a darkened room.

The **RUNNING** screen shows the feed rate, amount fed, flushing data if a feed & flush pump set is used, and the remaining volume to be delivered (remaining VTBD). The remaining VTBD is shown next to the  $\boxed{V}$  icon.

The **RUN** option can be selected from the following screens, if the proper parameters are programmed and the pump set is loaded:

- SET LOADED screens (Figures 7A-7D)
- ADJUST FEED screen, continuous mode (Figure 11A)
- ADJUST FLUSH screens (Figure 12A, 14A)
- HOLDING screens (Figure 18A)




Figure 15A. RUNNING mode screen for continuous feed mode.

Figure 15B. KTO.

#### Keep Tube Open (KTO) (Not Available in EZMODE)

The Keep Tube Open (**KTO**) function within the Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump will help prevent a blockage of the feeding set by moving the feeding solution periodically within the tube. On regular intervals, calculated from the pump's presently programmed feed rate, a small and clinically insignificant portion of feeding solution (0.2mL/interval) is moved through the tube. The amount of delivered fluid is monitored, tracked, and calculated into the final feed volume totalizer amounts displayed to the user. The **KTO** feature can be activated via the **RUNNING** screen. To activate the **KTO** feature, press the input button to the left of the **KTO** icon, see Figure 15A, to immediately input a desired time in minutes for which the pump will remain in the "Keep Tube Open" mode before resuming the previously programmed feeding regimen. Each press of the input button will increment the **KTO** time by 5-minute increments, up to 30 minutes, see Figure 15B. After "30" is reached, each additional button press will increment in larger time steps, up to a maximum of 240 minutes. Another button push after 240 will cycle the **KTO** back to 5 minutes. No other input is needed to start the countdown, as the time immediately counts down after the incremental value is reached.




Figure 16. Program Locked screen, which indicates the pump is in LOCK SETTINGS mode.




Figure 17. Padlock indicates pump is in LOCK SCREEN mode.

#### Lock Screen Feature (Not Available in EZMODE)

The Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump can be set to **RUN** and then locked out from unauthorized use. While the **LOCK SETTINGS** feature prevents only changes to the input parameters, the LOCK SCREEN feature will lock out any button presses, including RUN or **HOLD**, except power off.

Use **LOCK SCREEN** when the pump is running and there is a need to prevent tampering such as putting the pump into an unwanted **HOLD**. The **RUNNING** screen will show a padlock to indicate that no changes can be made without unlocking the screen.

To lock the screen, use the following procedure while the pump is in **RUNNING** mode:

• Press and hold the third > button from the top for at least 3 seconds, and then release it.

To unlock the screen, use the same procedure above.

## Hold

Pump operation can be paused by pressing > HOLD while it is *RUNNING*, such as in Figure 15A. The status line will indicate HOLDING and the yellow LED indicator will be lit up to give a quick visual indication of an informational status in a darkened room.

Note: the pump's alarm will sound if it has been left HOLDING for 10 minutes without further input.

There are five options available during *HOLDING*. See Figures 18A and 18B.




Figure 18A. *HOLDING* mode options, continuous feeding mode.

Figure 18B. *HOLDING* mode options, intermittent feeding mode.

#### **Clear Volume Fed**

Press > CLEAR VOLUME to re-zero the output feed totals "mL FED," "BOLUSES FED," and "mL FLUSH."

#### **Adjust Settings**

Press > **ADJUST SETTINGS** to return to one of the opening menus shown in Figures 7A-7D, thus allowing any setting to be modified. If the pump set has been removed after going into **HOLD**, the opening menu will show **LOAD A SET**, as seen in Figures 5A-5B.

#### **RESUME IN** \_\_\_\_ **MIN** (Not Available in EZMODE)

Press "**RESUME IN** \_\_\_\_**MIN**" once to immediately place the Kangaroo<sup>™</sup> ePump<sup>™</sup> into a 30-minute resume running mode. Pressing the **RESUME IN** \_\_\_\_**MIN** button again will cycle the "resume in" time back to zero. No other input is needed to start the countdown, as the time immediately counts down after the 30 minute value is reached.

#### **Run Mode**

Press > **RUN** to immediately restart pump operation. This option will not be available if the pump set is removed during **HOLD**.

#### **View History** (Not Available in EZMODE)

Press > **HISTORY** to access the *HISTORY* screen. Use this screen to review feeding and flushing history. See subsection "More Options" below for further information.

## Flush Now (Not Available in EZMODE)

During the *RUNNING* mode, the pump can be diverted to immediately perform a flush (flush on demand) by pressing > **FLUSH NOW** in the *RUNNING* menu, input a desired flush volume, then press > ENTER. See Figure 15A. When > FLUSH NOW is pressed, the screen will prompt for the volume to be immediately flushed, and will show the previously programmed flush volume as default but may be adjusted to the users desired volume. The screen will appear as the standard flush input screen as shown in Figure 12B (or 14B). Setting the **FLUSH NOW** flushing volume is temporary and will not change the main/periodic flushing volume previously programmed. If > **FLUSH NOW** is pressed by accident; the user may either press the cancel button or allow the screen to return to the **RUNNING** screen after 10 seconds, if there is no input for the flush volume.

If stopping a **FLUSH NOW** before the set volume has been delivered is desired, press **FLUSH NOW** and enter a volume of zero and press **FLUSH NOW** will stop immediately and the pump will return to the **RUNNING** mode.



Figure 19. MORE OPTIONS menu.

## **More Options**

Press > **MORE** in one of the opening menus, Figures 5A-5B or Figures 7A-7D, to access the **MORE OPTIONS** screen, Figure 19.

#### **Buzzer Alarm Volume**

Press > **BUZZER** in the **MORE OPTIONS** menu to access the **BUZZER** screen and increase or decrease the alarm buzzer volume. The alarm will be audible as the volume setting is changed. This sets the volume level for any warning or error that may occur with the pump. High Priority alarms, such as System Error, Rotor Error, and Pump Set Dislodged will always use the maximum volume setting and will be louder than other alarms.

#### **History** (Not Available in EZMODE)

Press > **HISTORY** in the *MORE OPTIONS* menu to access the *HISTORY* screen. Press the > buttons in Figure 20 to set the time history of interest. The totals of **FEED mL** and **FLUSH mL** will be shown for the requested history time.

**Warning:** A maximum of 72 hours of past history can be retrieved, excluding data from the current hour. **Note:** The 72-hour history includes time when the pump is powered OFF.

#### Language Selection

Press > **LANGUAGES** in the *MORE OPTIONS* menu to change the preferred language for the screen interface.

#### **Continuous/Intermittent Mode** (Not Available in EZMODE)

Press > **CONTINUOUS**/**INTERMITTENT** in the *MORE OPTIONS* menu to access the *SELECT MODE* menu, Figure 10B. Use this menu to select between continuous feeding mode and intermittent feeding mode. See subsection "Selecting Continuous or Intermittent Feed Mode" for more information.

Press > DONE in the MORE OPTIONS menu to exit the menu.




Figure 20. Feed and Flush totals for up to 72 hours of previous history, excluding data from the current hour.

## **BioTech Menu**

The **BioTech** screen will show technical information about the pump, and has a lock feature to lock the pump settings against unauthorized changes.

## Accessing the BioTech Screen

To get the **BioTech** screen, press the power button to start the pump, and press and hold the top ▶ button while the startup sequence (Kangaroo icon "hops" across the screen) executes. You must be pressing the top ▶ button at the end of the sequence in order to access the BioTech menu. Figure 21 shows the **BioTech** screen.

SN: <b>C</b> SU: FLASH: CERT. INTV: BEGION CODE:	12345678 1.02 1.05 6 1
BIOTECH OPTIO	ns B

Figure 21. BioTech screen.




Figure 22. BioTech option screen, with Lock Setting, EZMODE, and Communications features activated.

## **BioTech Screen Data**

SN is the serial number of the pump.

SW shows the software version number in this pump.

FLASH shows the version number of the flash data installed in this pump.

**CERT INTV** indicates the number of times the pump has been powered on since certification.

## **BioTech Options**


#### Lock Settings Feature — Program Locked

To "lock" the settings against unauthorized changes, press > LOCK SETTINGS to get the check mark, as shown in Figure 22. The settings can only be changed again by pressing > LOCK SETTINGS to remove the check mark. If the settings are locked and a change is attempted, the screen will indicate **PROGRAM LOCKED**, and will not allow access to the parameters, Figure 16.

The pump can be **PRIMED** or set to **RUN**, **HOLD**, etc. when settings are locked.

#### **EZMODE**

Press **EZMODE** to activate the **EZMODE** feeding option. A check mark will appear next to the **EZMODE** button to indicate that the feature is active. Pressing the **EZMODE** button again will remove the check mark indicating that the feature is off. If the **EZMODE** feature is turned on, the **EZMODE** icon will appear in the lower right corner of the LCD next to the battery icon (see Figure 9). See page 12 for **EZMODE** details.



## COMMUNICATION

Press > **COMMUNICATION** to activate the communication feature. A check mark will appear next to the > **COMMUNICATION** button to indicate that the feature is active. Pressing the > **COMMUNICATION** button again will remove the check mark indicating that the feature is off. If the **COMMUNICATION** feature is turned on, the **COMMUNICATION** icon will appear in the lower right corner of the LCD next to the battery icon (see Figure 23).

Turning the **COMMUNICATION** feature on will initiate data transfer of various pump settings, status, errors, etc. through the serial I/O port. If feature is desired, please contact technical support for additional information.

**Communication icon** 

Kangaroo™ ePump™ Enteral Feed and Flush Pump with Pole Clamp

Note: If EZMODE and the COMMUNICATION feature are both activated the icons for both features will alternate.

Press ► **DONE** to exit the **BioTech** menu. To re-enter this menu screen, power down the pump and then power it up again and hold the top ► button, as described above.

## **Section VI – Performance Evaluation**

#### System Performance Tests

A series of tests can be performed to verify pump performance. It is recommended that tests be run at least once every two years, or any time the pump is suspected of having improper performance. The test procedure and certificate are included on the User Manual CD as a separate document.

#### **Certifying the Accuracy of Pump Flow Rate**

Use the following procedure to check the flow rate accuracy:

- Perform test at room temperature 22° C  $\pm$  2° C (72° F  $\pm$  3° F) with a new pump set.
- Fill a new feed-only pump set bag with water to the 500 mL mark.
- Load the pump set.
- Suspend the pump set bag so that the bottom of the bag is 18 inches above the pump.
- Make sure pump is plugged into A/C power (do not run test with battery power).
- Program a Continuous Mode feed rate (see below for suggested example rates).
- Run the pump for 15 minutes to achieve steady state operation.
- Then begin collecting the water in a calibrated measuring container, such as a large-capacity graduated cylinder, for exactly 30 minutes.

**Note:** The amount of water collected in 30 minutes will be half of the value programmed for the hourly flow rate, within  $\pm$  10%, or 0.5 mL, whichever is larger.

Rate Set	Amount Expected in 30 Minutes
75 mL/hr	33.7 mL to 41.3 mL
150 mL/hr	67.5 mL to 82.5 mL

If the amount collected is outside the range, load a new pump set and re-run the test procedure to confirm results.

## Verifying functionality of the alarm system

The System Performance Test is sufficient to verify proper pump function. However, another quick test that can be performed to confirm audibility and function of alarms is:

1. Load a new feeding set onto the pump

#### 2. Leave the feeding set empty!

- 3. Run Auto Prime
- 4. Once the pump begins priming, it will issue a Feed or Flush Error alarm
- 5. Confirm that the audible alarm and colored LEDs all properly indicate a Feed or Flush Error

Caution: To allow proper operation of audible alarms, avoid blocking or obstructing the series of buzzer holes located in the back of the pump.

## Section VII — Cleaning

Cleaning should be performed as needed. It may also be desirable to define cleaning intervals based on knowledge of the environment in which the pump is used. Only personnel trained in the cleaning of medical devices should perform cleaning.

Caution: Do not immerse pump or A/C power adapter in water or other cleaning solution; clean using a damp (not wet) cloth or sponge. Failure to follow the cleaning procedures described herein could result in hazards to users. As with any AC powered electrical device, care must be taken to prevent liquid from entering the pump to avoid electrical shock hazard, fire hazard, or damage to electrical components.

If any of the following events occur, **do not use** the pump until it has been properly cleaned and serviced by personnel trained in servicing Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pumps:

- wetting of the pump's A/C power adapter or leakage into the pump interior during cleaning.
- spillage of large amounts of formula onto the pump exterior, or any spillage onto the A/C power adapter.

#### **General Cleaning Directions**

Cleaning of Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pumps must be performed as follows:

# Caution: Disconnect pump from AC power source before cleaning. After cleaning, do not connect to AC power source until pump and A/C power adapter are thoroughly dry.

- A mild detergent should be used for general cleaning. For disinfection a 10% hypochlorite solution, isopropyl alcohol or ammonium chloride based solutions may be used. Allow excess moisture to evaporate from the cord prior to use of A/C Power Adapter.
- **Do not use** strong cleaners such as Spray Nine<sup>™</sup>\*, pHisoHex<sup>™</sup>\*, Hibiclens<sup>™</sup>\*, or Vesta-Syde<sup>™</sup>\* or cleaners that may contain Ethoxylated C9-C11 alcohols, 2-Butoxyethanol, Hexachlorophene, Chlorhedidine Gluconate, Subtilisins, Sodium Tetraborate Decahydrate, Triethanolamine, or other similar ingredients because damage to the pump case housing can result.

## **Directions for Cleaning Pump Housing**

- Refer to General Cleaning Directions before starting.
- Clean outside surface with a damp (not wet) cloth or sponge and keep pump in upright position as much as possible.
- Avoid excess moisture near pole clamp area.
- Do not allow cleaning solution to enter the vertical vents in the housing, underneath the handle.

## **Directions for Cleaning Pump A/C Power Adapter**

- Refer to General Cleaning Directions before starting.
- Unless soiling is observed, the A/C power cord should not be cleaned.
- If cleaning of the power adapter is necessary, unplug from outlet and wipe the exterior surfaces of the wall plug with a cloth dampened with isopropyl alcohol.

#### Caution: Avoid exposing A/C power adapter to excess moisture, as this can lead to an electrical shock or fire hazard.

## **Directions for Cleaning Rotor Assembly**

- Refer to General Cleaning Directions before starting.
- Open the blue door enclosing the pump set loading area.
- Loosen rotor set screw with 5/64" (2 mm) allen wrench and gently pull rotor off shaft. After removing rotor, avoid getting any moisture in the rotor shaft opening.
- Clean rollers thoroughly with warm soapy water, or isopropyl alcohol if necessary.
- Be sure all parts of rotor are completely dry before putting it back onto shaft.
- To replace rotor, align set screw on hub of rotor with the flattened section of the output shaft.
- Push rotor into place and tighten set screw (Do not over tighten).

#### **Preventative Maintenance**

This pump must be periodically tested to assure proper functioning and safety. The recommended service interval is every 2 years. Performance testing may be done at the user's Biomedical Engineering Department, an outside service, or by Covidien Factory Service. To arrange for Covidien Factory Service in the US, call 1-800-962-9888. In Canada, call 1-877-664-8926. Outside of the US and Canada, please contact your local Customer Service for more information.

If a pump malfunctions, please contact your Covidien Representative or call Customer Service for instruction.

## Section VIII — Battery Replacement

See Figure 24 and do the following to replace the battery pack:

- Power down the unit.
- Loosen the battery door screw, on back of pump.
- Open battery door and disconnect wire harness. A small straight-bladed screw driver may be needed to disconnect the connector tab.
- Slide battery pack out of compartment.
- Slide a new battery pack into compartment with wire harness extending out.
- Align connector with red wire at upper right.
- Push connector in until properly latched.
- **Caution:** Ensure the battery wires are routed through the slot between the connector and the battery to prevent the wires from becoming pinched.
- Re-tighten the battery door screw; if there is any resistance while securing the door, recheck the wires to verify they are properly routed through the wire slot.




Figure 24. Installation of Battery Pack.

#### Notes:

- Replacement battery packs are available from Covidien (See Section XIII Service Part Numbers).
- The battery is packed partially charged. Charge battery 6 hours prior to use off A/C power.
- To assure the safe and reliable operation of this pumping system, do not attempt to operate without an approved battery pack installed in the pump.
- Dispose of used battery cells in accordance with local and facility guidelines.
- Keep away from children.

# Warning: Use of a battery pack not supplied by Covidien can be dangerous and voids all warranties and performance specifications.

The battery will charge continuously whenever the pump is plugged into a wall outlet. Eight hours of charging is required to fully recharge the battery pack.

For operational errors, the Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump display will show the error indicator icon, with specific information regarding the error condition, and will sound the alarm.

## The following conditions will trigger error conditions and will activate audible and visible alarms:

- Feed or flush containers are empty
- · Feed or flush tube becomes occluded between bag and pump
- Tube becomes occluded between pump and patient
- Battery low (alarm beeps continuously, about 15 minutes before shut-off)
- Unit is left in HOLD mode longer than 10 minutes without input
- Pump set tubing is improperly loaded around the rotor
- System errors, see below
- Feeding complete
- The MISTIC connector is removed during AUTOPRIMING or RUNNING

## **Priority Handling of Alarms**

The Kangaroo<sup>™</sup> ePump<sup>™</sup> has alarms that are broken into two different priorities: High Priority and Medium Priority. In all cases, High Priority alarms are the most important and override any other alarm conditions. System Error overrides other High Priority alarms. Other alarms are Medium Priority. These alarms all have the same priority except Low Battery alarm, which has higher priority than other Medium Priority alarms.

Below is the listing of alarm priorities for the pump:

#### **High Priority**

0: System Error Alarm

1: All Other High Priority Alarms

#### **Medium Priority**

2: Low Battery Alarm

3: All Other Medium Priority Alarms

## **LED Indicator Lights**

The pump status LED Indicator lights on the upper right of the pump provides a quick visual indication of the pump status, especially in darkened rooms.

A solid green light indicates the pump is ready for operation or feeding.

A solid yellow light indicates the pump is on hold or a medium priority alarm.

A solid or blinking red light indicates a high priority alarm.

## **Caregiver Alarm Notification**

All alarms are intended to be heard by operators that are within hearing range of the pump buzzer. The pump buzzer is located on the back of the pump. The pump is designed so the alarm can be heard within the patient's room, at a minimum. The display and LED alarm indicators are intended to be seen by an operator within the room, facing the front of the pump. Since audible alarms are limited by distance, it is recommended that the operator conduct a check to determine at what distance the alarm can still be heard.

**Note:** Going outside of the patient's room may make hearing the alarms more difficult.

The Error, Warning, and Information Screens are described below:

## System Error (High Priority Alarm)

The system error screen is the most general form of error, Figure 25. The only way to exit from a System Error is to power down. No **CONTINUE** option is allowed because of the severity of the error. An error number is also displayed on the screen, for reference purposes. This number should be reported when calling the Customer Service Support Line. Detection of this alarm condition may take up to 1 minute to occur under normal operating conditions.




Figure 25. General System Error screen. See Customer Service section.

Error #	Related Sub-system
0	Internal timer stopped.
1	Flash Chip missing OR version incompatible.
2	Too many graphic screens in memory.
4	Never got a ADC ISR for motor current.
5	Never got a ADC ISR for battery level.
6	Never got a ADC ISR for battery temperature.
7	Never got a ADC ISR for Ultrasonic Upstream Occlusion or Bag-Empty checks.
8	Never got a ADC ISR for Ultrasonic during Downstream Occlusion.
9	RS232 Communications Transmission Error.
10	Flash Write Error.
11	Flash Read Error.
12	General Stopcock Time Out.
13	Time Out for Buzzer Test
16	Flash Screen Read Error.
19	Flash Languages Read Error.
20	Could not find Stop Cock Flush Position during prime routine.

- 21 Could not find Stop Cock Feed Position during prime routine.
- 22 Motor Error during Tubing AutoLoad.
- 23 Motor Failsafe Time Out Error. Possible damaged gearbox encoder or flex cable.
- 25 Display Lockup Time Out Error.
- 33 Motor Control Failure.

Turn the power off and back on to try to clear the error. If the error cannot be resolved in a timely manner, press **U POWER** to stop operation of the pump and put a different pump into service.




Figure 26. Hold Error screen.

#### Hold Error (Medium Priority Alarm)

The **HOLD ERROR** screen will appear if the pump has been inactive, without input, for more than 10 minutes. See subsection "Hold" for description of *HOLDING* mode, Figure 26.

Press > **CONTINUE** to return to the previous screen. Where settings can be adjusted, the pump can be set to run immediately, or the pump can be set to run in a specified number of minutes. See Figures 18A-18B to see the options on the **HOLDING** menu screen.

If the error cannot be resolved, press **U POWER** to stop operation of the pump and put a different pump into service. Detection of this alarm condition may take up to 10 minutes to occur under normal operating conditions.

#### **Rotor Error** (High Priority Alarm)

The **ROTOR ERROR** screen appears during **RUNNING** or **PRIMING**, when the pump has detected that the pump set is no longer providing sufficient tension on the rotor. See Figure 27. The **ROTOR ERROR** generally results from a problem with the Pump Set tubing around the rotor.

Check that the pump set is not damaged, and re-load the tubing around the rotor as described in subsection "Loading Pump Sets". Press > **CONTINUE** to restart.

If the error cannot be resolved, press **U POWER** to stop operation of the pump and put a different pump into service. Detection of this alarm condition may take up to 70 minutes to occur at 1 mL/hr or up to 1 minute to occur at flow rates greater than 125 mL/hr.




Figure 27. Rotor Error screen. Re-load the set tubing and press CONTINUE.



Figure 28. Feed Error screen.

Feed Error (Medium Priority Alarm)

The **FEED ERROR** screen appears when the enteral formula is no longer being delivered because the bag is empty or because of a clog between the pump and the bag. The pump determines this by measuring the amount of fluid in the tube over time. Large amounts of foam or bubbles in the feeding solution can also be a cause for this alarm. Check the bag to see if it is empty and re-fill the bag as required. If the bag still contains feeding solution, remove cassette and check the bag side tubing for excessive foam or bubbles. Clear bubbles from line and reload the feeding set or replace with a new feeding set. Check the feed line to find the occlusion causing the blockage. If the occlusion cannot be cleared, load a new pump set, prime it, and restart the feeding, Figure 28.

If the error still cannot be resolved, press **O POWER** to stop operation of the pump and put a different pump into service. Detection of this alarm condition may take up to 360 minutes to occur at 1 mL/hr or up to 18 minutes to occur at flow rates greater than 125 mL/hr.

**Note:** An occlusion may pressurize the feeding set, which can result in an unintended bolus of feed or flush solution when the occlusion is cleared. This volume is less than 1ml.

### Flush Error (Medium Priority Alarm)

The **FLUSH ERROR** screen appears when the flushing solution is no longer being delivered because the flush bag is empty or because of a clog between the pump and flush bag. The pump determines this by measuring the amount of fluid in the tube over time. Large amounts of foam or bubbles in the feeding solution can also be a cause for this alarm. Check the bag to see if it is empty and re-fill the bag as required. If the bag still contains feeding solution, remove cassette and check the bag side tubing for excessive foam or bubbles. Clear bubbles from line and reload the feeding set or replace with a new feeding set. Check the flush line to find the occlusion causing the blockage. If the occlusion cannot be cleared, load a new pump set, prime it, and restart the feeding, Figure 29.



Figure 29. Flush Error screen.

If the error still cannot be resolved, press **U POWER** to stop operation of the pump and put a different pump into service. Detection of this alarm condition may take up to 2 minutes to occur.

**Note:** An occlusion may pressurize the feeding set, which can result in an unintended bolus of feed or flush solution when the occlusion is cleared. This volume is less than 1ml.




Figure 30. Flow Error screen.

#### Flow Error (Medium Priority Alarm)

The **FLOW ERROR** screen appears when the feeding solution is no longer being delivered because of a clog between the pump and the patient. The pump determines the presence of a clog by checking to see if fluid can be pumped away from the sensor below the feeding set valve while the valve is closed. Check the line to find the occlusion causing the blockage. If the error cannot be resolved, check the valve pocket in the pump set loading area for dampness or dirt, which may be a cause for a false error. Clean and dry the valve pocket. If the error still cannot be resolved, load a new Pump Set, prime it, and restart the feeding, Figure 30.

If the error still cannot be resolved, press **O POWER** to stop operation of the pump and put a different pump into service. Detection of this alarm condition may take up to 240 minutes to occur at 1 mL/hr or up to 14 minutes to occur at flow rates greater than 125 mL/hr.

**Note:** An occlusion may pressurize the feeding set, which can result in an unintended bolus of feed or flush solution when the occlusion is cleared. This volume is less than 1ml.

## Pump Set Dislodged Error (High Priority Alarm)

The **PUMP SET DISLODGED** screen will appear if the black ring retainer (MISTIC) is not properly loaded in the MISTIC pocket in the pump set loading area. Check and correct the positioning of the MISTIC retainer if possible. If the error cannot be resolved, load a new pump set, prime it, and restart the feeding, Figure 31.

If the error cannot be resolved, press **U POWER** to stop operation of the pump and put a different pump into service.



Figure 31. Pump Set Dislodged Error screen.




Figure 32. Battery Low error screen. Indicates that the battery must be recharged immediately.

## Figure 32.

**Battery Low** (Medium Priority Alarm)

Plug the A/C power adapter in to begin recharging. The pump will automatically return to the screen that was active prior to the error. The battery will charge continuously whenever the pump is plugged into a wall outlet. Eight hours of charging is required to fully recharge the battery pack.

The **BATTERY LOW** screen appears and the alarm beeps continuously when the battery needs to be recharged. There is approximately 30 minutes of battery life remaining when this screen appears,

#### Feeding Complete (Information)

The **FEEDING COMPLETE** information screen appears after completion of the programmed feeding, Figure 33. Press > **Power Down** to turn off the pump. Press > **CONTINUE** to return to the **SET LOADED** opening menu, Figures 7A-7D.



RUNNING				
SET USE	FEED RATE			
КТО	`125₩			
	423 mL			
	FED FLUSH <b>200</b> mL			
HOLD				
HOLD	EVERY 2 hrs			
FLUSH NOW	1000mL 🚆			

Figure 34. Running screen with the > 24 HRS of Pump Set Usage indicator in upper left of the screen.

LED Array

## Pump Set Use >24 Hours Warning

The pump set usage warning indicator will blink on the **RUNNING** screen if a pump set has been used for 24 or

more hours (hours actually running). It is recommended to replace pump sets after this length of usage. This icon is only an informational message and does not require action.

Figure 34 shows an example of the **RUNNING** screen, where a feed & flush set has been in use for more than 24 hours. The informational warning indicator shows up in the upper left on the screen.

Fix formatting. Low Battery Error is another item in the list. Bold it and add comma, no colon.

notification screen.

The pump status LED array on the upper right of the pump gives a quick visual indication of the pump status, especially in darkened rooms.

A green light status indicates normal pump operation.

A yellow light status indicates an informational or warning situation for a low battery condition, the pump in *HOLDING* mode, or one of the following medium priority alarms: LOW BATTERY ERROR: **HOLD ERROR**, **FEED ERROR**, **FLOW ERROR**, and **FLUSH ERROR**.

A red light status indicates kigh priority alarm. The screen will show one of the following error messages: SYSTEM ERROR, PUMP SET DISLODGED, or ROTOR ERROR.

Insert 'a'.

Kangaroo<sup>™</sup> ePump<sup>™</sup> Enteral Feed and Flush Pump with Pole Clamp

Table of Contents

## **Specifications**

#### **TUV Listing**



#### **Medical Electrical Equipment**

Kangaroo<sup>™</sup> ePump<sup>™</sup> Enteral Feed and Flush Pump with Pole Clamp

(1) Classified with respect to electrical shock, fire, and mechanical hazards in accordance with ES60601-1:2012, UL60601-1, IEC 60601-1:2012.

(2) Classified with respect to electrical shock, fire, mechanical and other specified hazards in accordance with CAN/CSA C22.2 No. 60601-1:08.

#### **Type Infusion Device** Volumetric

#### **Pumping Mechanism Rotary Peristaltic**

#### Pump Sets

Kangaroo<sup>™</sup> ePump<sup>™</sup> Feed-Only Set or Feed & Flush Set (with MISTIC Connector)

#### **Feeding Formula Delivery Rate**

1-400 mL/hr in 1 mL increments

## Feeding Formula VTBD

1-3000 mL in 1 mL increments

**Bolus Volume** 1-3000 mL in 1 mL increments

#### **Number of Boluses** 1-99

**Bolus Interval** 1-24 hours in 1-hour increments

#### **Flushing Solution Dose Range** 10-500 mL in 1 mL increments

Skip a line

**Flushing Solution Interval Range** 1-24 hr in 1 hr increments

#### Accuracy

 $\pm 10\%$  or 0.5 mL/hr, whichever is larger, with bag at 46 cm (18") above pump, at room temperature 22° C ( $\pm$  2° C (72° F  $\pm$  3° F), using water and a new pump set with less than the recommended 24 hours of maximum usage.

"Single Fault Condition Effect on Accuracy - In the event of a single fault short circuit in the pump electronics, a maximum bolus of 1.67 mL of extra fluid delivery may occur during feeding and 8.19 ml during flushing or priming."

#### **Feeding Set Occlusion Pressure**

15 psi (103 kPa) Nominal

#### Dimensions

Height: 16.8 cm (6.6") Width: 16.3 cm (6.4") Depth: 11.7 cm (4.6")

Kangaroo<sup>™</sup> ePump<sup>™</sup> Enteral Feed and Flush Pump with Pole Clamp

Remove quotes

Table of Contents

#### Weight

1.1 kg (2.4 lbs), 1.2 kg (2.7 lbs) with pole clamp

#### Material

Housing: ABS/Polycarbonate blend Pump set door: Polyester/ Polycarbonate blend

#### Power

Use AC power adapter for wall outlet usage. Pump operates on 9V DC, 1.5 Amps.

#### Battery

New, fully charged NiMH Battery Pack delivers  $\geq$  15 hrs at 100 mL/hr feed rate. Approximately 15 minutes prior to complete battery discharge, a low battery alarm will occur (see subsection "Battery Low" in Section IX). When complete discharge occurs, the pump will automatically turn itself off.

The battery will charge continuously whenever the pump is plugged into a wall outlet. Eight hours of charging is required to fully recharge the battery pack.

#### Alarms

- System Error
  Rotor Error
- Flow Error
  Feed Error
- Pump Set Dislodged Error
  Flush Error
- Hold Error
  Battery Low

#### **High Priority Alarm Volume**

Minimum of 64 dBA at 1 meter

#### **Medium Priority Alarm Volume**

Minimum of 47 dBA at 1 meter

Alarm Volume cannot be turned down below 47 dBA."

#### **Operating Temperature**

10° - 40° Č (50° - 104° F) 75% R.H. non-condensing

#### **Transport & Storage Temperature**


0° - 50° C (32° - 122° F) 95% R.H. non-condensing

Type of Protection Against Electrical Shock Class II, Internally-powered Equipment

**Degree of Protection Against Electrical Shock** Type BF

#### Mode of Operation Continuous or programmed intermittent operation

**Degree of Protection Against Ingress of Fluids** Drip-proof-IPX1



## Symbols on Pump and Accessories

The following symbols are found on the pump or on the accessories, such as the A/C adapter.

Symbol	Definition	Symbol	Definition	Symbol	Definition	
STERILE EO	Sterile using ethylene oxide	X	Dispose of as Electrical and Electronic Waste	Use- by	Use by or expiration date	
NON STERILE	Non-sterile	<b>SA</b> XX12345	UI recognized component certification symbol	LOT	Batch code	
<b>X</b>	This product does not contain natural rubber latex.		Caution: For indoor use only	SN	Serial Number	
Rx	Federal (USA) law restricts this device to sale by or on the order of a physician	<b>T</b>	Type BF applied part	M	Date of manufacture	
DEHP	DEHP Free		Class II Equipment (degree of protection against electrical shock), Double Insulated	$\bigcirc$	Power Symbol	
	Do not use if package is opened or damaged		Functional Ground	$\square$	Buttons for Selecting Screen Prompts	
i	Consult instructions for use	IPX1	Drip Proof (Degree of protection against ingress of fluids)		Power Source 9 VDC 1.5 A	
	Caution, consult accompanying documents	IPX0	Not water resistant	RoHS	RoHS	
MR	MR unsafe (magnetic resonance)	(((••)))	Non-ionizing electromagnetic radiation		DB9 Serial I/O Communication Port (Certification and Programming)	
$\diamond$	Flushing Solution	~	Alternating Current	D	UL Demko's Registered Safety Mark	
	Feed		Direct Current (DC)		Set Loading Diagram	
(2ah)	Do not use for greater than 24 hours		NRTL test mark; "Type Test" performed for compliance	(IV)	Level of Efficiency for External Power Supplies	
32°F 0°C	Operating Temperature	<b>CE</b> 0123	CE Mark	(ap		
32°F 0°C	Avoid temperature extremes	EC REP	Authorized representative in the European Community	White transition Stepped Connector	Assure the white ENFit transition stepped connector is securely tightened. If applicable, assure the cap is securely tightened.	
0% 95%	Humidity	REF	Catalog number	Enteral Feed Catheter \		
<b>Ť</b>	Keep Dry	<b>~~~</b>	Manufacturer			

The circuitry of the Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump is not customer serviceable. In particular, electronic assembly rework by non-authorized Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump technicians will likely affect accuracy.

All service personnel must be properly trained and qualified with operation of the Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump. Improper service may impair operation of the pump.

#### **Return for Repair**

Call Customer Service for an Authorized Return Number and shipping instructions, using the appropriate phone number below:

**United States** Covidien Phone: 1-800-962-9888 **Canada** Covidien Phone: 1-877-664-8926

For all customers outside of the U.S. and Canada, please contact your local customer service.

## Section XII — Maintenance

For general maintenance issues not discussed below, contact Customer Service (Section XI).

**Warning:** Do not open the main housing, as there are no user-serviceable parts inside. Opening of device may affect function of device and voids the warranty.

**Warning:** Do not perform maintenance when the pump is in operation.

The following maintenance items/parts can be replaced by the customer on the Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump. See Section XIII for Service Part Numbers and Section XI for contacting Customer Service.




Figure 35. Blue Door Replacement.

## Side Door enclosing Pump Set Loading area

- To replace the side door, see Figure 35 and do the following:
- 1. Place top door hinge onto the top hinge pin.
- 2. Carefully press down on the bottom of the door to spread the door hinges apart.
- 3. Place the bottom door hinge onto the bottom hinge pin.

## **Battery Pack**

See Section VIII for full instructions on replacing battery pack. During long-term non-usage or storage (greater than 9 months) of the pump, battery discharge and deactivation may tend to occur. It is recommended that the battery be disconnected during this time. The battery may need to be charged and discharged several times to restore optimal capacity.

## **Battery Door**

To replace battery door, see Figure 36 and do the following:

- 1. Align the battery door hinges with the long hinge pins.
- 2. Snap the hinges onto the hinge pins.
- 3. Screw down the door to hold in place.

## **A/C Power Adapter**

See Section IV for initial setup, including A/C Power Adapter attachment.




Figure 36. Battery Door Replacement.


## **Pole Clamp**

See Section IV for initial setup, including attachment of the pole clamp to the pump.

## Section XIII – Service Part Numbers

To place an order for parts, or if technical assistance is required, call customer service.

The Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump contains a limited number of serviceable parts, Figure 37. User maintenance is to be performed only by appropriately qualified technical personnel.






Figure 37. User Serviceable Parts.

## North/South America

772055 Kangaroo™ ePump™ 500 mL Pump Set 773656 Kangaroo™ ePump™ 1000 mL Pump Set 674668 Kangaroo™ ePump™ 500 mL Feed Set with 500 mL Flush Bag 773662 Kangaroo™ ePump™ 1000 mL Feed Set with 1000 mL Flush Bag 716154 Kangaroo™ ePump™ Sterile 100 mL Burette Set 775659 Kangaroo™ ePump™ Safety Screw Spike Set 775100 Kangaroo™ ePump™ Safety Screw Spike Set with 1000 mL Flush Bag 776150 Kangaroo™ ePump™ 100 mL Burette Recertification

## UK

7771067 Kangaroo™ ePump™ 1000ml Pump Set 7770647 Kangaroo™ ePump™ 500ml Pump Set 7781047 Kangaroo™ ePump™ 1000 mL Feed Set with 1000 mL Flush Bag 7781617 Kangaroo™ ePump™ 1600 mL Feed Set with 1000 mL Flush Bag 7774017 Kangaroo™ ePump™ ENPlus Spike Set 7774027 Kangaroo™ ePump™ Dual ENPlus Spike Set 7774037 Kangaroo™ ePump™ ENPlus 3 in 1 Pump Set 7774057 Kangaroo™ ePump™ ENPlus Spike Set with 1000 mL Flush Bag 7774067 Kangaroo™ ePump™ ENPlus 3 in 1 1000 mL Feed Set with 1000 mL Flush Bag

#### EMEA

777401 Kangaroo™ ePump™ ENPlus Spike Set 777402 Kangaroo™ ePump™ Dual ENPlus Spike Set 777403 Kangaroo™ ePump™ ENPlus 3 in 1 Pump Set 777405 Kangaroo™ ePump™ ENPlus Spike Feed Set with 1000 mL Flush Bag 777406 Kangaroo™ ePump™ ENPlus 3 in 1f Feed Set with 1000 mL Flush Bag 777407 Kangaroo™ ePump™ ENPlus Spike feed with vented Spike flush 777408 Kangaroo™ ePump™ ENPlus 3 in 1 Feed Set with Flush Bag

#### **Limited Warranty:**

- Covidien warrants to the original purchaser ("Customer") that this newly manufactured enteral feeding pump ("Pump" or "Pumps") will be free of defects in materials and workmanship, under normal use, for three (3) years from the date of shipment from Covidien. This Limited Warranty as applied to pump batteries and power cords is limited to one (1) year from the date of shipment from Covidien for all pumps.
- 2. This Limited Warranty does not extend to routine maintenance of pumps such as cleaning and all recommended performance tests set forth in this pump operation and service manual which remain the sole responsibility of Customer. Failure of Customer to perform cleaning, routine maintenance and recommended performance testing on any pump as outlined in this pump operation and service manual may void this Limited Warranty.
- 3. Customer agrees that, with the exception of customer serviceable parts and troubleshooting steps outlined in this pump operation and service manual, Covidien or its authorized dealer must perform pump repairs.
- 4. This Limited Warranty does not cover any pump, product or part that:
  - (a) has been operated in an unsuitable environment or used for purposes other than intended;
  - (b) has been subjected to unauthorized or non-Covidien repair or use of non-Covidien supplied parts;
  - (c) has been altered, misused, abused or neglected;
  - (d) has been subjected to fire, casualty or accident;
  - (e) suffers damage caused by Customer's negligent acts or omissions; or
  - (f) suffers damage beyond normal wear and tear.
- 5. For purposes of this Limited Warranty, "damage beyond normal wear and tear" includes without limitation:
  - (a) Damage to housing, LCD, display overlay or power supply;
  - (b) PCBA damage due to fluid ingress;
  - (c) Use of non-qualified power supply or battery; or
  - (d) Use of unauthorized cleaning fluids.
- 6. If a pump does not operate as warranted during the applicable warranty period, Covidien may, at its option and expense,
  - (a) repair or replace the defective part or pump; or,
  - (b) refund to Customer the purchase price for the defective part or pump.
- 7. Dated proof of original purchase is required to process warranty claims. Removal, defacement or alteration of serial lot number voids this Limited Warranty.
- 8. Shipping costs for pumps being returned to Covidien shall be borne by Customer. Customer is responsible for proper packaging for return shipment. Loss or damage in return shipment to Covidien shall be at Customer's risk.
- 9. Covidien disclaims all other warranties, expressed or implied, including any implied warranty of merchantability or fitness for a particular purpose or application other than as expressly set forth in the product labeling. In no event shall Covidien be liable for any incidental, indirect or consequential damages in conjunction with the purchase or use of the Pump, even if advised of the possibility of the same.

The Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump with pole clamp has been built and tested according to UL 60601-1 (2nd Edition), ES60601-1:2012, CAN/CSA C22.2 No. 60601-1-08, and EN60601-1-2 Standards.

The Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump with pole clamp is intended for use in the electromagnetic environment specified below. The user of the pump should assure that it is used in such an environment.

The Kangaroo™ ePump™ Enteral Feed and Flush Pump with Pole Clamp is intended for use in the electromagnetic environment specified below. The user of the Kangaroo™ ePump™ Enteral Feed and Flush Pump with Pole Clamp should assure that it is used in such an environment.					
<b>Emissions Test</b>	Compliance	Electromagnetic Environment - Guidance			
RF emissions (CISPR 11)	Group 1	The Kangaroo™ ePump™ enteral feed and flush pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.			
RF emissions (CISPR 11)	Class B	The Kangaroo <sup>™</sup> ePump <sup>™</sup> enteral feed and flush pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.			
Radiated Disturbance Immunity (EN60601-1-2 / IEC 61000-4-3: 2006 + A1: 2007 + A2: 2010)	Complies				
Conducted Disturbance Immunity (EN60601-1-2 / IEC 61000-4-6:2013)	Complies				
Power Frequency Magnetic Field Immunity (EN60601-1-2 / IEC 61000-4-8:2009)	Complies				
Voltage dips and sags Immunity (EN60601-1-2 / IEC 61000-4-11:2004)	Complies				
Electrical Fast Transient / Bursts Immunity (EN60601-1-2 / IEC 61000-4-4:2012)	Complies				
Electrostatic Discharge Immunity (EN60601-1-2 / IEC 61000-4-2:2008)	Complies				
Surge Immunity (EN60601-1-2 / IEC 61000-4-5:2005 +Cor 1: 2009)	Complies				

Guidance and Manufacturer's Declaration – Electromagnetic Immunity	
Surdanice and manufacturer 5 beenaration Electroniughette initiality	

The Kangaroo™ ePump™ enteral feed and flush pump with pole clamp are intended for use in the electromagnetic environment specified below. The user of the Kangaroo™ ePump™ enteral feed and flush pump with pole clamp should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) (EN 61000-4-2 per EN 60601-1-2: 2015)	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	$\pm 2$ kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode	± 1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Nominal Mains Voltage (VNOM) 100 VAC and 240 VAC @ 50 Hz Voltage Dips 100% of VNOM for 10 mSec (0.5 Line Cycles) at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 100% of VNOM for 20 mSec (1 Line Cycle) at 0° 70% of VNOM for 500 mSec (25 Line Cycles) at 0° Interruptions 100% of VNOM° for 5000 mSec (250 Line Cycles)	No degradation of performance or loss of function.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Kangaroo™ ePump™ enteral feed and flush pump with pole clamp require continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field (EN 61000-4-8 per EN 60601-1-2: 2015)	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

#### Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The Kangaroo™ ePump™ enteral feed and flush pump with pole clamp are intended for use in the electromagnetic environment specified below. The customer or the user of the Kangaroo™ ePump™ enteral feed and flush pump with pole clamp should assure that it is used in such an environment

Immunity Test	IEC 60601 Test Level		Compliance Level	Electromagnetic Environment -Guidance	
Conducted RF IEC 61000-4-6			6 Vrms	6 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Kangaroo™ ePump™ entera
Radiated RF (EN 61000-4-3 per EN 60601-1-2: 2015)	Band        380-390 MHz        430-470 MHz        704-787 MHz        800-900 MHz        1700-1990 MHz        2400-2570 MHz        5100-5800 MHz	Í.	Modulation        PM, 18 Hz, 50% 27 V/m        PM, 18 Hz, 50% 28 V/m        PM, 217 Hz, 50% 9 V/m        PM, 217 Hz, 50% 28 V/m	27 V/m 28 V/m 9 V/m 28 V/m 28 V/m 28 V/m 9 V/m	any part of the Kangaroo <sup>™</sup> ePump <sup>™</sup> enteral feed and flush pump with pole clamp, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance Not applicable $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump is used exceeds the applicable RF compliance level above, the Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Kangaroo<sup>™</sup> ePump<sup>™</sup> enteral feed and flush pump.

<sup>b</sup>Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

# Recommended separation distances between portable and mobile RF communications equipment and the Kangaroo™ ePump™ Enteral Feed and Flush Pump with pole clamp

The Kangaroo<sup>TM</sup> ePump<sup>TM</sup> enteral feed and flush pump with pole clamp is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Kangaroo<sup>TM</sup> ePump<sup>TM</sup> enteral feed and flush pump with pole clamp can help prevent electromagnetic interference by maintaining the minimum distance between portable and mobile RF communications equipment (transmitters) and the Kangaroo<sup>TM</sup> ePump<sup>TM</sup> enteral feed and flush pump with pole clamp can help prevent electromagnetic interference for mobile RF communications equipment (transmitters) and the Kangaroo<sup>TM</sup> ePump<sup>TM</sup> enteral feed and flush pump with pole clamp can help prevent electromagnetic interference for mobile RF communications equipment (transmitters) and the Kangaroo<sup>TM</sup> ePump<sup>TM</sup> enteral feed and flush pump with pole clamp recommended below, according to the maximum output power of the communication equipment

Rated maximum output	Separation distance according to frequency of transmitter m			
power of transmitter				
W	150 kHz to 80 MHz d = 1,2√P	80 MHz to 800 Mhz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz d = 2,3√P	
0,01	0,12	0,12	0,23	
0,1	0,38	0,38	0,73	
1	1,2	1,2	2,3	
10	3,8	3,8	7,3	
100	12	12	23	

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.





Manual No. HP112480

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.

\*\*\*Spray Nine is a trademark of Spray Nine Corporation

™\*pHisoHex is a trademark of The Mentholatum Co.

<sup>™\*</sup>Hibiclens is a trademark of Regent Medical Ltd.

™\*Vesta-Syde is a trademark of Steris Inc.

Other brands are trademarks of a Covidien company.

© 2012 Covidien. Covidien IIc, 15 Hampshire Street, Mansfield, MA 02048 USA. ECREP Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.

REV 03/2020