Photo-Therapy 4000

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Contents

For Your Safety and that of Your Patients	17
General WARNINGS and CAUTIONS	18
Intended Use	18
System overview	19
Symbols	19
Preparation	20
Mounting on stand	20
Connecting to mains power supply	20
Using the unit	21
Operation	22
Safety precautions	22
Switching on the lamp	23
Care of patient	23
Monitoring service life	23
Shut-down	24
Cleaning	25
Stripping down	25
Disinfecting and cleaning	25
Fault – Cause – Remedy	26
Maintenance	26
Maintenance intervals	26
Replacing fuse	26
Technical Data	27
Disposal	28
Design and Operating Principles	29
Order List	29

Definitions

WARNING

A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or in damage to the equipment or other property.

NOTE

A NOTE provides additional information intended to avoid inconvenience during operation.

For Your Safety and that of Your Patients

WARNING

For a full understanding of the performance characteristics of this device, the user should carefully read this manual before use of the device.

Strictly follow these Instructions for Use

WARNING

Strictly follow these Instructions for Use. Any use of the product requires full understanding and strict observation of all portions of these instructions. The device is only to be used for the purpose specified under "Intended Use" on page 18 and in conjunction with appropriate patient monitoring. Observe all **WARNING** and **CAUTION** statements throughout this manual and all statements on device labels.

Maintenance

WARNING

The device must be inspected and serviced regularly by trained service personnel. Repair of the device may also only be carried out by trained service personnel. Dräger recommends that a service contract be obtained with DrägerService and that all repairs also be carried out by them. Dräger recommends that only authentic Dräger repair parts be used for maintenance. Otherwise the correct functioning of the device may be compromised. See chapter "Maintenance".

Accessories

WARNING

Only the accessories indicated on the order list have been tested and approved to be used with the device. Accordingly it is strongly recommended that only these accessories be used in conjunction with the specific device. Otherwise the correct functioning of the device may be compromised.

Not for use in areas of explosion hazard

WARNING

This apparatus is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

Safe connection with other electrical equipment

WARNING

Electrical connections to equipment which is not listed in these Instructions for Use should only be made following consultation with the respective manufacturers.

Patient safety

The design of the equipment, the accompanying literature, and the labeling on the equipment take into consideration that the purchase and use of the equipment are restricted to trained professionals, and that certain inherent characteristics of the equipment are known to the trained operator. Instructions, warnings, and caution statements are limited, therefore, largely to the specifics of the Dräger design. This publication excludes references to various hazards which are obvious to a medical professional and operator of this equipment, to the consequences of product misuse, and to potentially adverse effects in patients with abnormal conditions. Product modification or misuse can be dangerous.

Patient monitoring

The operators of the system must recognize their responsibility for choosing appropriate safety monitoring that supplies adequate information on equipment performance and patient condition. Patient safety may be achieved through a wide variety of different means ranging from electronic surveillance of equipment performance and patient condition to simple, direct observation of clinical signs. The responsibility for the selection of the best level of patient monitoring lies solely with the equipment operator.

Dräger Medical GmbH

General WARNINGS and CAUTIONS

The following WARNINGS and CAUTIONS apply to general operation of the device. WARNINGS and CAUTIONS specific to subsystems or particular features appear with those topics in later sections of the manual.

Note on EMC/ESD risk for the device function

General information on electromagnetic compatibility (EMC) pursuant to international EMC standard IEC 60601-1-2: 2001. Electromedical devices are subject to special precautionary measures concerning electromagnetic compatibility (EMC) and must be installed and put into operation in accordance with the EMC information included in the technical documentation which is available from DrägerService on request. Portable and mobile RF communications equipment can affect medical electrical equipment.

Intended Use

The Phototherapy Unit is designed to reduce the serum bilirubin concentration of premature babies and neonates by means of phototherapeutic radiation.


WARNING


Connector pins with an ESD warning sign should not be touched and no connections should be made between these connectors without implementing ESD protective measures. Such precautionary procedures may include antistatic clothing and shoes, the touch of a ground stud before and during connecting the pins or the use of electrically isolating and antistatic gloves. All staff involved in the above shall receive instruction in these ESD precautionary procedures.

System overview

- 1 Switch for the white lighting – area which is illuminated. A green indicator light in the switch also comes on.
- 2 Switch for the therapy – switches on the blue phototherapy lamps. A green indicator light in the switch also comes on.
- 3 Hour meter – records operating time for blue phototherapy lamps.
- 4 Handles for attaching the phototherapy lamp to the stand and for carrying, are situated to the left and right of the unit.


- 5 Mains inlet with built-in fuse drawer.
- 6 Mains switch selector (only for 110 to 115 V and 120 to 127 V phototherapy lamps).
- 7 Rating plate.

Stand

- 8 Arms to hold the phototherapy lamp. During normal use the lamp is secure when the handles are fitted into the arm notches.
- 9 Screw for adjusting the angle of the phototherapy lamp up to 30°.
- 10 Clamp to adjust the height of the phototherapy lamp between 1070 and 1620 mm (from the underside of the lamp).
- 11 Stand with castors for positioning the phototherapy lamp over an incubator or a warming cot.
- 12 Castors with brakes (where applicable) for locking the stand in position.

Symbols

Symbol Explanation


	Observe the Instructions for Use
	Observe minimum clearance of 30 cm
	Use light-tight eye protection
	ESD warning label
	Disposal information

Preparation

Mounting on stand

Mount the phototherapy lamp in the arms on the stand by the handles:

- 1 Insert into the front oval notches first
 - then slide in at an angle and
- 2 insert into the rear slanted notches.
 - Check that the unit is firmly in place.


Connecting to mains power supply

- 3 Mains voltage must correspond with operating voltage given on the rating plate.
 - The white switch for the white light and the blue switch for therapy on the front of the device must be in the »OFF« position.
- 4 Connect the mains cable to the phototherapy lamp and then to the mains.

For 110 to 115 V and 120 to 127 V phototherapy lamps

- 5 Check the position of the voltage selector switch. Choose the switch setting according to line voltage supplied, either 110 to 115 V or 120 to 127 V.

Using the unit

When used on the stand


WARNING

Minimum clearance between the lower edge of the phototherapy lamp and the patient = 30 cm. Risk of patient overheating.

Phototherapy in incubator

Example: Incubator 8000 NC


- 1 Loosen the clamp and adjust the height of the stand for the intensity of light required (see table, page 28).
 - Re-tighten the clamp.
 - Push the base of the stand under the incubator from the rear until the lamp housing is centred over the incubator.
- 2 Lock the castors (where applicable). The trolley can no longer be moved.

Phototherapy in warming cot with radiant heater

Example: Babytherm 8000 OC

- 3 Loosen the clamp and adjust the height of the stand for the intensity of light required (see table, page 28).
 - Re-tighten the clamp.
 - Push the base of the stand under the warming cot from the rear.
- 4 Turn the setscrew and tilt the Phototherapy Unit (max. 30° possible) until the light shines on the side of the patient. Radiant heat from the heater must reach the patient unimpeded.
- 5 Lock the castors (where applicable), so that the trolley can no longer be moved.

When used without the stand

The lamp may also be placed on the canopy of the incubator itself or on the canopy of the warming cot.

- Check the load bearing capacity of the canopy! Check the stability of the phototherapy lamp!
- Minimum clearance to the patient = 30 cm!
- Secure the lamp to prevent it from falling!

Operation

Safety precautions

CAUTION

Regularly measure the bilirubin levels of the patient.

WARNING

Cover the patient's eyes with light-tight eye protection to prevent damage to the eyes.

Beware of the possible toxic effects of photoisomers of bilirubin.

WARNING

Phototherapy can lead to unnoticed increased water losses. Ensure sufficient replacement, otherwise there is a risk of injury to the patient.

The duration of the therapy is determined by clinical requirements. Monitor the patient constantly.

If a radiation distance of less than 40 cm is required for therapeutic reasons, patient monitoring must be intensified.

Changes in ambient conditions, e.g. a draft, may influence the temperature balance of the patient.

WARNING

Additional heat sources, e.g. radiant heaters, increase the supply of heat to the patient and must be taken into account. Avoid exposing the patient to direct sunlight. Risk of patient overheating.

WARNING

The use of reflective foils can induce dangerous body temperatures.

WARNING

The absorption of light in the patient's skin leads to heat production. The core temperature of the infant can increase.

Use in the incubator:

Reduce the air temperature setting by about 2 °C and reduce the humidity setting.

Use of a radiant heater in manual mode:

Adjust the heat output of the radiant heater accordingly.

Strictly follow the specific Instructions for Use of the heater.

Monitor the patient's core temperature.

NOTE

Do not block the air vents of the unit.

WARNING

Whilst the unit is in operation, direct eye contact with the radiation source over a long period must be avoided. The eyes can be damaged through the effect of the radiation.

WARNING

Protect patients in the immediate vicinity of the phototherapy unit from the radiation produced, e.g. protective goggles should be worn and screens used. The eyes can be damaged through the effect of the radiation. Risk of sunburn.

WARNING


Users should avoid remaining in the radiation area for extended periods. Risk of sunburn.

WARNING

Whilst the unit is in operation, do not use any inflammable detergents or medicaments. Do not place medicaments or infusion solutions in the radiation zone. Fire risk.

Switching on the lamp

- 1 Push the white switch to switch the white light on (where applicable). The green indicator light in the switch will come on.
- 2 Push the blue switch to switch the blue phototherapy lamps on. The green indicator light in the switch will come on. It takes about 1 minute to reach maximum output.
- 3 The hour meter records the operating time of the blue phototherapy lamps.

If therapy is interrupted and the lights are turned off

- Phototherapy lamps must be allowed to cool down for about 1 minute before turning them back on.

Care of patient

- Monitor the infant's core temperature.
- Increase fluid supply to meet infant's increased fluid consumption, e.g. increase fluid supply by parenteral infusion.
- Check that eye protection remains in place at all times.
- Check serum bilirubin concentration periodically.

Monitoring service life

The intensity of the effective light decreases with use. After 500 operating hours the intensity is more than 85 % and after 1000 operating hours the intensity is more than 80 % of the initial value.

Hour meter shows total operating time for the phototherapy lamp.

Blue phototherapy lamps must be replaced every 1000 operating hours by trained service personnel.

WARNING

Use only the phototherapy lamps indicated on the Order List. The use of other radiation sources can impair the reliability and effectiveness of the phototherapy.

Shut-down

When therapy is finished:

- Switch off the blue phototherapy lamps.
- Switch off the white lights (where applicable).
- Disconnect from the mains power supply.
- Remove the phototherapy unit from the incubator/warming cot canopy.

When the mobile stand has been used:

- Release the castor brakes (where applicable).
- Move the stand away from the incubator/warming cot.
- Disconnect from the mains power supply.
- Re-adjust the angle of the lamp to the »zero« position.
- Loosen the clamp and lower stand to transport height.
- Store the phototherapy unit in a safe place.

Cleaning

NOTE

Even accessories designed to be reused (e.g. after cleaning) have a limited service life. Due to a number of factors connected with handling and cleaning (e.g. disinfectant residues can attack the material more intensely during autoclaving), increased wear can occur and the service life can be markedly shortened. If there are any external signs of wear such as cracks, deformations, discolorations, peeling materials or the like, these accessories have to be replaced.

Stripping down

- Lift the phototherapy lamp off the stand.
- Wipe off any obvious stains with a disposable cloth.

Disinfecting and cleaning

Use only preparations classified as "surface disinfectants" for disinfecting.

For material compatibility, we recommend preparations based on:

- aldehydes
- quaternary ammonium compounds.

Due to the chemical composition and direct effect on material compatibility, preparations based on the following components are of only limited suitability:

- halogen-releasing compounds
- strong organic acids
- oxygen-releasing compounds.

In general, the manufacturer's recommendations must be followed when selecting the preparation. The manufacturer is legally liable for the information concerning the range of applications and possible material damage caused by these preparations.

NOTE

Do not use any disinfectants or cleaning agents containing alcohol.

For users in the Federal Republic of Germany, we recommend the use of disinfectants listed in the current DGHM list (DGHM = German Society for Hygiene and Microbiology). The DGHM list (published by mhp-Verlag GmbH, Wiesbaden) also specifies the active basis of each disinfectant.

For countries where the DGHM list is unavailable, we recommend products based on the above active bases.

The following surface disinfectants can be recommended:

Dismozon [®] pur	Bode Chemie GmbH & Co, Hamburg
Incidur [®]	Henkel Hygiene GmbH, Düsseldorf
Sekusept [®] pulver	Henkel Hygiene GmbH, Düsseldorf

Follow the manufacturer's instructions for use.

Wipe-disinfecting the phototherapy lamp

- Wipe off obvious stains using a disposable cloth and detergent.
- Wipe-disinfect surfaces.
- After disinfecting, wipe surfaces again with a clean damp cloth and then dry thoroughly.
- Do not allow any liquid to get into the lamp.

CAUTION

Use only the recommended cleaning agents and disinfectants, as there is a risk of tension cracks in the acrylic glass if, for instance, alcohol is used.

CAUTION

Do not allow any liquids to penetrate the inside of the phototherapy unit or to come in contact with the light bulbs. Otherwise the correct functioning of the device may be compromised.

After each care/cleaning/disinfecting operation, check that the unit is functioning correctly.
(see page 23, "Switching on the lamp")

Fault – Cause – Remedy

Fault	Cause	Remedy
No lights come on	Not connected to mains power	Connect to the mains
	Faulty fuse	Replace fuse
	Not switched on	Switch on
Some pairs of lights do not come on	Faulty switching unit	Have switching unit replaced by DrägerService or trained service personnel
	Faulty lights	Have lights unit replaced only by DrägerService or trained service personnel

Maintenance

Maintenance intervals

Clean and disinfect device or device parts before each maintenance step – and also when returning for repair!

Inspect device	annually by specially trained personnel
Clean vents	annually by specially trained personnel
Replace blue phototherapy lamps	every 1000 operating hours by DrägerService or trained service personnel

Replacing fuse


WARNING

Before replacing fuse, disconnect lamp from mains power supply. Danger of electric shock.

CAUTION

Fuses must be of the same value as those they replace, see "Technical Data" on page 27.

- Disconnect from the mains power supply.
- Disconnect the mains cable from the phototherapy lamp.
- 1 Open the fuse drawer below the mains inlet and remove the fuses.
- Replace the faulty fuse with a new one of the same value.
- Close the fuse drawer and re-connect to the mains.
- Check that the lamps work.

Technical Data

Ambient conditions

Data on temperature apply for stable conditions.

Unit temperature = ambient temperature

During operation:

temperature	15 °C to 35 °C (59 °F to 95 °F)
atmospheric pressure	700 to 1060 hPa
humidity	20 to 80 % relative humidity (without condensation)

During storage:

temperature	-10 °C to 60 °C (-4 °F to 140 °F)
atmospheric pressure	500 to 1060 hPa
humidity	0 to 80 % relative humidity (without condensation)

Operating data

	2M 21 000	2M 21 700
Mains voltage	220 to 240 V AC	110 to 115 V or 120 to 127 V (switchable)
Current consumption	0.6 A	1.2 A
Frequency	50/60 Hz	50/60 Hz
Power consumption	120 W	120 W
Fuses	2 x T 1.0 A; 250 V IEC 127-2/V	2 x T 1.6 A; 250 V IEC 127-2/V
Fluorescent lamps		
therapy (blue light)	4 x 18 W	4 x 18 W
lighting (white light)	2 x 18 W	2 x 18 W
Noise emission	none	none
Weight		
lamp	5.5 kg (12.1 lbs)	7.2 kg (15.9 lbs)
stand	13.6 kg (30 lbs)	13.6 kg (30 lbs)
Dimensions in mm		
lamp	470 x 285 x 130 mm (18.5 x 11.2 x 5.1 inches) (L x W x H)	
base of stand	690 x 600 x 100 mm (27.2 x 23.6 x 4 inches) (L x W x H)	
height adjustment of stand	1070 mm to 1620 mm (42 to 64 inches)	

Standards

	conforms to EN 60601-1
Electromagnetic compatibility	tested to EN 60601-1-2: 2001
Enclosure protection class	I
Classification as per EC Directive 93/42/EEC Annex IX	Class II a
UMDNS code	13-037
Universal Medical Device Nomenclature System – Nomenclature for medical products	

Light intensity on optical axis

The effective light intensity for bilirubin dissociation E_{bi} depends on how far away the lamp housing (lower edge) is. If a higher intensity of light is required, the white light lamps may be replaced with blue phototherapy lamps (by DrägerService or trained service personnel). The hours used by the two additional phototherapy lamps will **not** be registered by the hour meter.

Maximum values

Distance	cm		30	40	50
E_{bi}	mW/cm ²	4 blue light lamps	2.4 ±0.4	1.6 ±0.3	1.2 ±0.2
		6 blue light lamps	3.4 ±0.6	2.3 ±0.4	1.6 ±0.3

Distribution of light intensity

This measurement was performed at a vertical distance of 40 cm on a bed surface measuring 50 x 30 cm. Maximum value = 100 % Intensity specified as a percentage of the maximum value. Precision: ±5 %.

63	82	91	82	63	10 cm
67	88	100	88	67	10 cm
63	82	91	82	63	10 cm

10 cm 10 cm 10 cm 10 cm 10 cm

Average values on working surface given:

4 blue light lamps: 1.3 mW/cm²
 6 blue light lamps: 1.8 mW/cm²

Disposal

For countries subject to the EU-Directive 2002/96/EC:

This device is subject to EU Directive 2002/96/EC (WEEE). It is not registered for use in private households, and may not be disposed of at municipal collection points for waste electrical and electronic equipment. Dräger Medical has authorized a firm to dispose of this device in the proper manner: for more detailed information, please contact your local Dräger Medical organization.

Design and Operating Principles


The phototherapy lamp consists of a plastic-coated aluminium housing with electronic switching units, fluorescent lamps for light therapy and lighting, control circuits, fuses and an "hours used" meter.

The lamp may also be mounted on a mobile stand (order separately) which allows it to be used at an angled position when used with infra-red heated warming cots, such as the Dräger Babytherm 8000.

The light is helpful for malfunctions of bilirubin metabolism. If a newborn baby develops icterus because detoxification of the liver has been inadequate, therapy to lower bilirubin concentration may be required as early as the first day. Phototherapy has been found to be a reliable treatment for excessive bilirubin concentration.

The light converts the bilirubin into a water-soluble form by a photochemical reaction and so enables natural excretion via bile and urine.

The blue portion of the light spectrum is most effective in decomposing bilirubin at an absorption peak of 460 nm. The Dräger fluorescent lamps emit light in this peak.

Order List

Name and description	Order No.
Photo-Therapy 4000 220 to 240 V, with mains cable	2M 21 000
Photo-Therapy 4000 110 to 115 V and 120 to 127 V, switchable, with mains cable	2M 21 700
Stand for Photo-Therapy 4000	2M 21 190
Spare parts	
Fluorescent lamp for lighting (white)	2M 21 009
Fluorescent lamp for therapy (blue)	2M 21 010

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These Instructions for Use only apply to

Photo-Therapy 4000

with the Serial No.:

If no Serial No. has been filled in by Dräger, these Instructions for Use are provided for general information only and are not intended for use with any specific machine or device.

This document is provided for customer information only, and will not be updated or exchanged without customer request.

Directive 93/42/EEC
concerning Medical Devices

IM24

Manufactured by:

Dräger Medical GmbH

Moislinger Allee 53 – 55

D-23542 Lübeck

Germany

+49 451 8 82- 0

FAX +49 451 8 82- 20 80

<http://www.draeger.com>

Only for Canada:

Distributed by

Draeger Medical Canada Inc.

120 East Beaver Creek Road Suite 104

Richmond Hill Ontario L4B 4V1

+1 866 343 2273 (Toll Free)

+1 905 530 2387

FAX +1 905 763 1890

90 29 150 – GA 6141.421 de/en/fr/es/nl/sv/ru/it/ro/bg

© Dräger Medical GmbH

Edition: 9– 2012-12

Dräger reserves the right to make modifications to the equipment without prior notice.