

Baxter Sigma Spectrum IV INFUSION PUMP



The Baxter Sigma Spectrum IV Infusion Pump System (non-wireless model), for purchase or rental, was designed with features to help prevent pump-related medication errors when programming an IV infusion. SIGMA Spectrum is an easy-to-use pump with Dose Error Prevention Software and standard IV administration sets.

- PC-based, pharmacy edited and controlled, customized in-house library of all IV and epidural drugs, along with their safe delivery parameters.
- Configurable Dose Error Reduction Software (DERS) to help reduce pump-related adverse drug events
- Large and small volume parenterals (LVP).
- Uses standard gravity IV sets from (Hospira and Baxter)
- Removable, rechargeable battery pack for easy replacements
- Lithium Ion, 1800 mA/h, 7.4 VDC nominal. 12 hr. recharge time

*Product photo may include various optional features. Please consult with your sales associate for more details.

Call your local US Med-Equip representative or 877-677-7767 today for more information.

usmedequip.com