

## Respironics BiPap Synchrony/Synchrony S



**General Category:**

Bi-level Ventilation

**Manufacturer:**

Respironics (now owned by Phillips)

**Model Number:**

BiPap Synchrony/Synchrony S

**Customers Who Order:**

Homecare & LTAC

**Application of Equipment:**

The BiPAP Synchrony provides more natural and comfortable noninvasive ventilation specifically for COPD patients and others suffering from respiratory insufficiency.

**Specifications:**

Unit may have all or some of the following features: Continuous positive airway pressure, Timed- and Pressure-Control Modes. Digital Auto-Trak Sensitivity has enhanced leak compensation to allow for better matching of the patient's breathing pattern. Pressure range of 4 to 30 cm H<sub>2</sub>O and the ability to collect patient therapy, compliance and machine data using an optional built-in modem.

**Associated Disposables:**

Reusable Filter  
Ultra-Fine Filter  
Breathing Circuit  
Face Mask – Full (S)  
Face Mask – Full (M)  
Face Mask – Full (L)  
Face Mask – Nasal (S)  
Face Mask – Nasal (M)  
Face Mask – Nasal (L)