

US Med-Equip Stainless Steel Six-Leg Heavy Duty IV Pole Model 8060



Heavy Duty and Extra Stable Base

The 8060 heavy duty stainless steel IV pole was designed to support the largest IV pump on the market. Its sturdy construction also makes it suitable to accommodate all pole mounted equipment (pump not included).

Product Features

- ◆ Telescoping stainless steel pole adjusts from 73" to 99.5"
- ◆ Six-leg, sturdy cast iron base
- ◆ Removable pop-off roller blade casters
- ◆ Easy to attach urine bag hooks (optional)

Specifications

◆ Base Material	Cast Iron	◆ Max Height Inches	99.5 in
◆ Disposable	No	◆ Min Height Inches	73 in
◆ Foldable	No	◆ Number of Hooks	4, 8
◆ HPIS Classification	770_130_0_0	◆ Number of Legs	6
◆ Instrument Family	IV Pole	◆ Product Type	IV Stand, Caster, O2 Tank Holder, Rake Hook Top, Guide Handle
◆ Instrument Finish	Stainless Steel	◆ Sterile	No
◆ Latex Free	Yes	◆ Type of Hook	Rams Horn
◆ Material	Stainless Steel, Rubber, Metal, Chrome Plated Steel	◆ UNSPSC Code	42222104, 52141605, 42192208
◆ Urine Bag Hooks (optional)		◆ Unit Compatibility	Mds80600
◆ Replacement Removable Casters (optional)			

7028 Gessner Rd. ◆ Houston ◆ Texas, 77040 ◆ U.S.A.
877-677-7767 (toll-free) ◆ FAX: 713-983-8861

www.usmedequip.com



Equipped for success.