

GE Marquette MAC 5500 EKG Machine



Product Specifications

Manufacturer GE Medical Systems

Features:

- 12-/15-Lead High-End EKG System for Resting-EKG's and Stress-Test
- Easy-to-use applications and features streamline productivity and workflow
- Advanced algorithms set the standard in EKG analysis and interpretation
- Improving Clinical Trial Workflow with CT Data Guard
- Multiple Network Capabilities (Ethernet and wireless)
- Internal Modem
- Barcode and magnetic card scanning options
- XML data export capabilities
- Internal Memory for 200 EKG's, SD-card removable Media
- Easy updating via Optioncodes
- Seamless connectivity to the MUSE® cardiology information system speeds data storage and EKG retrieval

The MAC 5500 represents a new definition of EkG speed and accuracy!

Built on GE's proven innovation in EKG acquisition and analysis, the **MAC® 5500** is GE's premier EKG system, delivering advanced disease management capabilities through its industry-leading collection of standards-setting algorithms and advanced networking.

The MAC 5500 system offers the sophistication required for advanced EKG applications, while its ease of use extends this level of performance to the broadest range of users possible. And, it's part of the complete GE suite of networked, non-invasive testing solutions designed to maximize efficiency and productivity in hospitals, clinics, office based practices and research institutes.

*Product photo is from the manufacturer and may include various optional features. Please consult with your sales associate for more details.

9777 West Gulf Bank Rd., Suite #20 ♦ Houston ♦ Texas, 77040 ♦ U.S.A.
877-677-7767 (toll-free) ♦ FAX: 713-983-8861

www.usmedequip.com


Equipped for success.