

Chattanooga OptiFlex 3 Knee CPM



Product Specifications

Manufacturer	Chattanooga
Knee Flexion	120°
Knee Extension	-10°
Speed	30° - 150°/minute
Pauses	0 - 30 seconds
Force Reversal	35 - 75 lb (16 - 34 kg)
Maximum Patient Weight	350 lb (159 kg)
Calf Length Range (knee joint to sole of foot)	10" - 23.5" (25.4 - 59.7 cm)
Thigh Length Range	12" - 19" (30.5 - 48.3 cm)
Mains Power	100 - 240 V, 50/60 Hz
Unit Weight	27 lb (12 kg)

Leading Therapy In A New Direction!

The sleek design of the OptiFlex 3 CPM is therapeutically functional as well as appealing to the eye. The ultra-wide carriage is designed to accommodate the typical knee replacement patient, athletes and pediatric patients. Designed to provide anatomical motion for CPM Knee Patients from Pediatric to Adult.

Product Features

- Ultra wide carriage - 11" (28 cm) from center of carriage
- Independent delay time
- Patient lockout feature
- Universal Left / Right Design
- Easy-to-read backlit display
- Patient compliance/usage meter
- Quick conversion from adult to pediatric

*Product photo is from the manufacturer and may include various optional features. Please consult with your sales associate for more details.

9777 West Gulf Bank Rd., Suite #20 ♦ Houston ♦ Texas, 77040 ♦ U.S.A.
877-677-7767 (toll-free) ♦ FAX: 713-983-8861

www.usmedequip.com


USMED-EQUIP

Equipped for success.