



GE Corometrics 250CX Maternal / Fetal Monitor



(STOCK PHOTO FROM MANUFACTURER)

Comprehensive perinatal monitoring that enables you to efficiently care for both mother and baby.

The Corometrics* 250cx Series maternal/fetal monitoring system is a full antepartum/intrapartum fetal monitor that includes innovative features that helps you deliver a high quality birthing experience and enhanced patient care. It monitors uterine and fetal activity including fetal heart rate, and has enhanced maternal monitoring features including temperature, oxygen saturation, ECG and non-invasive blood pressure. GE Healthcare offers two models within the Corometrics 250cx Series: the 259cx, which is the fully equipped maternal/fetal monitoring system, and the 256cx, which is a fetal monitor system.

The 250cx Product Features:

- **GE-exclusive Smart BP*** delays blood pressure measurement during contractions to reduce patient distraction and provide more meaningful diagnostic data.
- **Provides the light you need to see clearly**
 - ChartLight with backlit keys and a dimmable display helps you work comfortably while your patients remain relaxed in a low-light environment.
- **Allows for exceptional pick-up and minimal repositioning**
 - Our watertight transducers have 9-crystal ultrasound technology that allows for easier pick-up and less repositioning, even for difficult to monitor patients. The 9-crystal design achieves a wider focal region for more uniform coverage at greater depths than traditional 7-crystals.
- **Minimizes monitor interactions**
 - Since a majority of your monitor interactions can be done using our one-touch hot keys, you'll spend less time with the monitor and more time with patients.

Learn more about the technology that makes it easier to monitor vitals, such as:

- **Maternal blood pressure**
 - DINAMAP* SuperSTAT* non-invasive blood pressure technology speeds up determinations by as much as 66% after the first reading to enhance patient comfort and allow you faster access to patient data.
- **Maternal temperature**
 - Integrated Exergen Temporal Arterial Thermometer allows you to obtain maternal temperature using non-invasive infrared technology. Readings automatically print on the fetal strip.
- **Maternal oxygen levels**
 - The monitor provides a choice of integrated Masimo® SET or Nellcor® maternal pulse oximetry, making it easy for you to monitor, document, and even view the waveform.
- **MECG**
 - Selectable 3 Lead maternal ECG allows clinicians to provides clinicians QRS wave amplitude and meaningful diagnostic data.
- **Additional technology includes:**
 - **Centricity* Perinatal (CPN)** compatibility, which streamlines the capture and archiving of patient data and helps speed workflow and ensure data availability.
 - **Telemetry Interface** allows for continuous monitoring without changing transducers during ambulation or laboring in water.