

GE Carescape V100 Vital Signs Monitor



Product Specifications

Manufacturer	GE Carescape
Height	7.7 in (19.5 cm)
Width	8.6 in (21.9 cm) without temperature 10.0 in (25.4 cm) with temperature
Depth	5.3 in (13.5 cm)
Weight	5.4 lb (2.4 kg) including battery
Mountings	Self-supporting on rubber feet or pole mounted
Portability	Carried by recessed handle

On the go mobility. Clinical excellence and the quality you expect!

The CARESCAPE V100 monitor includes the same advanced parameters and algorithms as other higher acuity GE monitors, so you can consistently measure across all care areas:

- GE DINAMAP SuperSTAT, Classic or Auscultatory blood pressure algorithms support speed, comfort and artifact rejection
- Choose from three SpO2 measurements including GE Ohmeda, Nellcor or Masimo
- Three different temperature options to choose from: 1) Exergen TemporalScanner, a non-invasive, fast, easy and accurate thermometer, does not come into contact with mucous membranes, 2) Alaris Turbo Temp provides fast oral or rectal temperatures at the bedside, as well as axillary temperatures in monitor mode in as fast as 7 seconds. And 3) Alaris Tri-Site thermometer provides fast oral, rectal and axillary temperatures at the bedside. It provides predictive measurements with greater tolerance to probe placement in the vicinity of the sublingual pocket in as fast as 11 seconds.

*Product photo is from the manufacturer and may include various optional features. Please consult with your sales associate for more details.

7028 Gessner Rd. ♦ Houston ♦ Texas, 77040 ♦ U.S.A.
877-677-7767 (toll-free) ♦ FAX: 713-983-8861

www.usmedequip.com

USMED-EQUIP®
Equipped for success.