

Philips M2601A Telemetry Transmitter



The following Philips transmitters can be used with the Philips Telemetry System:

- . standard - ECG and SpO2
- . standard - ECG only
- . EASI - ECG and SpO2
- . EASI - ECG only

To aid in identification, standard ECG transmitters have dark green labels and EASI transmitters have purple labels.

The Philips Transmitter (EASI and standard ECG version) is battery powered and worn by the patient. It acquires the patient's ECG and SpO2 signals (if available), processes them, and sends them via the antenna system to the receiver. Measurements are then displayed at the Philips Information Center. The transmitter can also be connected via an infrared link to the Wave Viewer or TeleMon to provide display of patient measurements and waveforms at the patient's side.