

## Pulmonetic Systems LTV 1200 Portable Ventilator



### Product Features

#### ◆ Anytime transport:

Turbine technology in the ventilator eliminates the need for a high-pressure air source, allowing transport ventilation where and when you need it.

#### ◆ Optimizable trial settings:

By optimizing trial settings for each of your patients, you can ensure favorable levels of support throughout the weaning process.

#### ◆ Oxygen resource management tool:

The O2 Conserve feature reduces oxygen consumption as necessary, extending the sometimes scarce resource of oxygen during transport or emergency situations. It may be turned off when oxygen resources return to full supply.

### ***Excellent performance and safety features makes this the ideal portable ventilator!***

The LTV® Series 1200 ventilator supports infant, pediatric and adult patients in ICU, ER and patient transport with invasive or noninvasive ventilation pre-sets. These settings can be easily refined using the touch-turn-touch interface on the LED display. The ventilator also provides a wide range of ventilation therapies to meet demanding patient needs, including volume control, pressure control, pressure support and spontaneous breath types. Combined with the spontaneous breathing trial function, the ultra-sensitive flow trigger facilitates weaning.

\*Product photo is from the manufacturer and may include various optional features. Please consult with your sales associate for more details.

9777 West Gulf Bank Rd., Suite #20 ♦ Houston ♦ Texas, 77040 ♦ U.S.A.

877-677-7767 (toll-free) ♦ FAX: 713-983-8861

[www.usmedequip.com](http://www.usmedequip.com)

  
\*\*\*\*\*  
Equipped for success.