

Gynnie®

OB/GYN Stretcher



Adjustable footrests feature simple, one-handed operation

4" thick Ultra Comfort mattress provides unparalleled support

Sealed mattress and clean lines prevent fluid collection

Compact design with pop-up push handles and retractable fifth wheel improves maneuverability

Intuitive Patient Positioning

Gynnie combines the trademark durability and features of a standard Stryker stretcher with supreme OB/GYN functionality. Ideal for pelvic exams and general use, this multi-functional stretcher can outfit an entire ER department.



Glide footrests out, position patient and lower foot section to provide excellent pelvic access.



Gynnie can be configured with optional OB/GYN accessories.



Medical

3800 E. Centre Ave.
Portage, MI 49002
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Standard Features

- 4" (10 cm) Ultra Comfort Mattress with sealed seams
- 8" (20 cm) Omni Surface™ casters with wheel covers
- 26" (66 cm) wide patient surface
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics with constant rate descent
- Dual side-mounted foot control hydraulics
- Retractable fifth-wheel steering
- Four-wheel, steel-ring brakes with dual side and head activators
- Glideaway™ siderails
- Integrated storage tray with oxygen bottle holder
- Permanently-attached folding foot section
- Permanently-attached foot supports
- Pop-up push handles
- Roller bumpers
- Three IV receptacles
- Unique one-handed foot section release lever

Specifications

| | | |
|---------------------------------------|-------------|---------------|
| • Model Number | 1061 | |
| • Overall Length | 81.5" | (207 cm) |
| • Overall Width | | |
| Siderails Up | 31.5" | (80 cm) |
| Siderails Down | 26.75" | (68 cm) |
| • Weight Capacity | 500 lbs | (228 kg) |
| • Height Range (to litter top) | | |
| High | 36" | (91 cm) |
| Low | 22.5" | (57 cm) |
| • Litter Positioning | | |
| Backrest | 0° - 90° | |
| Foot Section | 0° - 90° | |
| Trend./Reverse Trend. | ±18° | |
| • Patient Surface | 26" x 77" | (66 x 196 cm) |
| • Siderails | 13" x 39" | (34 x 99 cm) |
| • Caster Diameter | 8" | (20 cm) |

Stryker reserves the right to change specifications without notice.

Warranty

- One year parts, labor and travel or two years parts

Optional Features

- C-spine cassette film holder
- Calf supports and brackets
- Defibrillator tray/footboard/extender
- Backrest X-ray cassette holder
- Hoop for drainage basin
- IV caddy
- IV poles
 - Permanent folding two- or three-stage
 - Tethered
 - Standard removable
- Patient restraints
- Serving/instrument tray
- Siderail pads
- Upright oxygen bottle holder
- Extended warranty available