

## CareFusion Alaris Signature 7230 Dual Channel Infusion Pump



### Product Specifications

Manufacturer	CareFusion
Overall Dimensions	10.5" (w) x 8.6" (h)
Weight	8.4 lbs.
Alarms	Air in Line, Battery Depleted, Channel Malfunction, Computer Link Failure, Flow Sensor Unplugged, Hold Time Exceeded, Instrument Malfunction, Key stuck, Latch Open, No Upstream Flow Detected, Occlusion Downstream, Occlusion Upstream, Primary Flow Detected during
Power Supply	100- 120 V~,50/60 Hz

### ***Excellent performance and safety features makes this the ideal infusion pump!***

The Alaris Signature 7230 is dual-channel infusion pump with independent programming on each channel. It can be used for general IV infusions, antibiotics, blood infusions, enteral infusions, TPN and hyperalimentation, and lipids. Monitoring options provide increased awareness of occlusion detection and battery charge condition. It features a dual-rate piggyback function and panel lock for improved safety.

A multi-step program function allows a sequential program to deliver up to nine rates and volumes; can also be programmed using volume and time. Multiple infusions can be preprogrammed over a period of up to 24 hours; delay start and loading dose features.

\*Product photo is from the manufacturer and may include various optional features. Please consult with your sales associate for more details.

9777 West Gulf Bank Rd., Suite #20 ♦ Houston ♦ Texas, 77040 ♦ U.S.A.  
877-677-7767 (toll-free) ♦ FAX: 713-983-8861

[www.usmedequip.com](http://www.usmedequip.com)

\*\*\*\*\*  
Equipped for success.