


NEW!

smiths medical  
bringing technology to life

# Medfusion® 4000 Syringe Infusion Pump with PharmGuard® Infusion Management Software Suite

Now with  
wireless capability




*Infusing Confidence  
with every medication delivered*

medfusion

TM

## Confidence delivered by...

- ✓ Improved pump functionality
- ✓ Improved library
- ✓ Wireless or wired data collection and report generation



### Closing the CQI Loop for Syringe Infusion

The PharmGuard® infusion management software suite promotes continuous quality improvement by allowing customized drug library configurations to be created, stored, then deployed to the Medfusion® 4000 pumps.



### Medfusion® 4000 Syringe Infusion Pump with PharmGuard® Medication Safety Software

- Built on the same Medfusion® pump platform known and trusted by clinicians
- +/- 2% accurate delivery at 0.01 mL/Hr to 1130mL/Hr flow rates
- Compatible syringe list is updateable
- Extended memory for history data storage
- Wireless / wired connectivity for device updates and CQI reporting
- IPX 3 fluid resistance

### PharmGuard® Toolbox Medication Safety PC Software

- Drug library configurations created for easy deployment to pumps
- Navigation improvements based on user feedback
- Configurable clinical alerts and Quick Library support best practices at bedside
- Allows easy migration from previous configurations
- Easy management of compatible IV and enteral syringes
- Windows®-based software databases for easy data reporting and data access

### PharmGuard® Server Software

- Deploys and monitors drug library configuration and firmware updates
- Provides Clinical / IT / Biomed reports based on collected CQI data
- Supports wired and wireless (802.11b/g) communications
- Uses industry standard OS, database, systems and protocols
- Flexible configuration options to meet hospitals' stringent IT needs

**Smiths Medical ASD, Inc.**  
St. Paul, MN 55112, USA  
Phone: 1-214-618-0218  
Toll-Free USA: 1-800-258-5361  
[www.smiths-medical.com](http://www.smiths-medical.com)

Find your local  
contact information at:  
[www.smiths-medical.com/  
customer-support](http://www.smiths-medical.com/customer-support)

Smiths Medical is part of the global technology business Smiths Group plc. Product(s) described may not be licensed or available for sale in all countries. Medfusion, PharmGuard and the Smiths Medical and Medfusion design marks are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. Rx only. ©2011 Smiths Medical. All rights reserved. 10/11 IN193825

**smiths medical**