

Cadd Prizm 6101 Ambulatory PCS Pump



General Category:
Ambulatory Pump/PCS

Manufacturer:
Cadd Prizm

Model Number:
6101

Customers Who Order:
Homecare & Home IV Companies

Application of Equipment:
The Cadd Prizm Multi-therapy Ambulatory Pump is a compact epidural pain management device used for transport administration of medication.

Specifications: Units may have all or some of the following features and modes: 4 delivery modes: PCA, Continuous, TPN, and Intermittent

Optional Accessories:

Bolus Cable	(Pt. Unknown)
Keys	(Pt. 21-2303-24)
Plastic Carry Case	(Pt. 21-2487)

Associated Disposables:

IV Set
Cassette

Provided By Customer
Provided By Customer